

STOP PRESS ** STOP PRESS **** STOP PRESS **** STOP PRESS **** STOP PRESS**

Company: Helapet Ltd
Contact: Ruth Hewitt - 01582 501983
Date: 14th April 2009
Title: New Berner non-latex gloves for handling cytotoxics

The long awaited non-latex version of the Berner Glove for handling cytotoxics is here at last!

The new Berner neoprene gloves are soft and comfortable to wear combining the lengthy breakthrough times associated with manufactured textiles together with the stretchy qualities needed to fit onto an isolator cuff.

This superb combination will give in excess of 480 minutes breakthrough time for Cisplatin and 90 minutes for carmustine at the same time offering excellent grip and tactile sensitivity.

The gloves are pliable enough to stretch over the cuff of an isolator and are strong to resist tearing, making them an ideal choice for aseptic units that are latex free.

The gloves are manufactured in compliance with EN 374-3:2003 and in line with PPE directive 89/686/EEC.

To try these exciting new gloves in your department, please call **0800 0328 428** for more information and **FREE SAMPLES**

