

PRESS RELEASE ** PRESS RELEASE **** PRESS RELEASE**

Title: New Chill Pouch added to Helapet Porter range

From: Helapet Ltd
Contact: Ruth Hewitt – Product Manager

Date: 25th March 2008

New Chill Pouch added to Helapet Porter Range

A tamper evident Chill Pouch is now available to transport temperature sensitive drugs from pharmacy manufacture to end user.

The new A4 Chill Pouches are fully validated and will comfortably accommodate 5 full 50ml syringes, maintaining the temperature between 2 and 8 deg C for in excess of 8 hours.

The Chill Pouches are easy to use and are supplied complete with medicool pack, divider and absorbent pad. The pouches are secured using the sealing tape and can easily be individually identified by using a document wallet.

For more information, please contact our sales team on 0800 0328 428 or visit our website www.helapet.co.uk



Ends