

HELAPET

'Filter Devices, Medical & Cleanroom Disposables'

To: The Editor
Date: 19 January 2002
Subject: Helapet 'STERICLEAN WIPE'
From: Glenn Kenworthy,
Marketing Manager, Helapet. Ltd.
Tel: 01296 632244 e.mail: glennk@helapet.co.uk

For Immediate Release

PRESS RELEASE

Helapet STERICLEAN WIPE – The world.wide.wipe.

Helapet introduces a **large individual sterilised alcohol wipe** in response to requests received from customers. The new Helapet **STERICLEAN WIPE** is ideal for cleaning large surfaces in your clean room/clean cabinet or isolator. It is individually packed in a convenient sachet, allowing for ease of use and precise stock control.

“Many of our customer’s have been looking for a large **sterilised** alcohol wipe for some time,” state’s Glenn Kenworthy, Marketing Manager at Helapet. “We were unable to source a product that really satisfied their requirements, so we decided to create our own. We are sure the Helapet **STERICLEAN WIPE** is the best large alcohol wipe available on the market today”

Each **STERICLEAN WIPE** is (approx.) 20 x 22cm and impregnated with 70% Isopropyl Alcohol BP before being sterilised by gamma irradiation. The product is packaged as 10 packs of 25 pieces (250 pieces) per box, allowing pharmacy, laboratory and cleanroom staff the ability to easily manage their usage.

The Helapet sales team can provide full product information and support, free samples and further information on 01296 632244 or you can e.mail us at: info@helapet.co.uk

www.helapet.co.uk

ends.