



# Helapet

## Pharmacy Workmat

Constructed of low-linting and absorbent cellulose material with an impermeable polyethylene backing, Helapet Pharmacy Workmats deliver essential drip and spill protection during drug preparation activities. Available in three convenient sizes, Helapet Pharmacy Workmats can be used in confined areas including cleanroom isolators and prep-trays.

Prepared in-house within an ISO Class 7 cleanroom environment, Helapet supply this product individually poly-bagged and gamma sterile for smooth aseptic transfer. A non sterile larger size is also available.



### Features and Benefits

- ✓ Highly absorbent materials for drips and spills
- ✓ Low linting materials
- ✓ Range of sizes to suit specific applications
- ✓ Packed in a heat sealed polybag
- ✓ Sterile and non sterile options

### Ordering Information

#### WP4657

46 x 57cm Helapet Pharmacy Workmat non sterile • 50 pcs/pck

#### XX2025

20 x 25cm Helapet Pharmacy Workmat sterile • 20 x 10 pcs/pck

#### XX2846

28 x 46cm Helapet Pharmacy Workmat sterile • 10 x 10 pcs/pck



Triple bagged packaging available

Enquire today for a quotation

Enquire today for  
**FREE**  
samples



## Technical Data

Material	Laminated cellulose	
Dimensions (w x l)	WP4657	46 x 57cm
	XX2025	20 x 25cm
	XX2846	28 x 46cm
Basis Weight (g/m)	159	
Thickness @ 53kPA (µm)	290	
Absorbency	0.75L/m <sup>2</sup>	
Packaging	WP4657	50 x pieces per pack
	XX2025	Individually poly bagged*; 10 x polybags per sealed outer bag, 20 x outer bags per pack
	XX2846	Individually poly bagged*; 10 x polybags per sealed outer bag, 10 x outer bags per pack
Compatibility	ISO Class 4 and Grade A environments	
Sterilisation	Gamma irradiation	
Shelf Life	Non Sterile	5 years from date of manufacture
	Sterile	3 years from date of manufacture

\*triple bag packaging available, enquire for a quotation

