



20mm Male Luer Lock (MLL) Vial Adapter

Helapet 20mm MLL Vial Adapter provides a secure connection between female luer lock (FLL) devices and a standard 20mm diameter vial.

This product constitutes part of the Vial2Bag® IV Transfer System*, allowing for direct vial reconstitution with any standard IV bag.

Using a refined plastic spike for easy insertion, the shrouded adapter design protects the user from the risk of accidental injury.

Presented individually sterile in moulded packaging, the MLL Vial Adapter can be easily fastened using non-touch technique, reducing the risk of cross contamination.

*Vial2Bag® transfer system (part ref: IV70101) available from Helapet.



IV72024

Features and Benefits

- ✓ Safe, stable application to 20mm vials
- ✓ Compatible with Vial2Bag® Transfer System
- ✓ Refined plastic spike for needle-free access
- ✓ Individually blister packed
- ✓ Sterile
- ✓ CE marked

Ordering Information

IV72024

20mm Male Luer Lock (MLL) Vial Adapter
sterile • 150 pcs/pck



Enquire
today for
FREE
samples



Technical Data

Material	Polycarbonate
Diameter	20mm
Blister	Polyethylene Terephthalate Glycol (PETG) Copolymer
Luer access	Accepts ISO 594 Compliant Female Luer Lock (FLL) devices and Halkey Roberts® needlefree connections
Packaging	Individually blister packed; 150 x blisters per pack
Sterilisation	Gamma irradiation
Compliance	CE marked; 510(k) cleared by United States Food and Drug Administration
Shelf Life	5 years from date of manufacture

