

A4 Cleanroom Paper



Helapet A4 Cleanroom Paper is made from quality cellulose and is coated to provide a smooth surface to minimise release of particles and fibres in cleanroom environments. Helapet Cleanroom Paper can retain toner from copiers and laser printers to reduce smudging and is compatible with both aqueous and oil-based ink.

Features and Benefits

- ✓ Smooth surface paper to minimise contamination
- ✓ Retains toner from copiers and laser printers to reduce smudging
- ✓ Compatible with aqueous and oil-based ink
- ✓ Double bagged for cleanliness



MP0003

Ordering Information

MP0003

A4 • non sterile • 2,000 sheets/pck

Technical Data

Dimensions	210mm x 297mm (UK A4 Size)
Material	Coated cellulose
Colour	Blue
Weight per sheet	105g
Packaged	10 x double bagged reems of 200 sheets
Sterility	Non sterile

Enquire today for
FREE
samples

Helapet Limited | Lister House | Blackburn Road | Houghton Regis | Bedfordshire | LU5 5BQ

Web: helapet.co.uk | **Phone:** +44 (0)1582 501980 | **Email:** sales@helapet.co.uk

Helapet Limited is a ClniMed Group Company.