



# Ampoule Breaker

## Sterile and non sterile



A tough and versatile polypropylene container, the Helapet Ampoule Breaker facilitates the safe breaking of ampoules, greatly reducing the risk of injury. A lightweight and portable device, the opening slot can easily accommodate different ampoule neck sizes.

With the ability to contain up to 200 ampoule caps, the Ampoule Breaker is ideal for mobile medical services and areas of limited waste disposal. Users can securely fasten the opening with an integrated mouth closure, ready for convenient disposal.

Available non sterile, or double bagged and gamma sterile for smooth aseptic transfer into high grade cleanroom environments.



AB7000

### Features and Benefits

- ✓ Small and compact design
- ✓ Safely holds up to 200 caps
- ✓ Integrated mouth closure between uses
- ✓ Instructions for use imprinted on the front
- ✓ Sterile and non sterile




XX7000

### Ordering Information

**AB7000**  
Ampoule Breaker  
non sterile • 48 pcs/pck

**XX7000**  
Ampoule Breaker  
double bagged • sterile • 25 pcs/pck



Triple bagged  
packaging available



Enquire  
today for  
**FREE**  
samples



## Clinical Waste and Spillage

### Technical Data

Material	Polypropylene	
Dimensions (mm)	H130 x L70 x D40	
Weight (g)	60.6	
Colour	Yellow	
Properties	Integrated mouth closure and printed instructions for use	
Packaging	Non sterile	48 x units loose packed per pack
	Sterile	Individually double poly bagged*; 25 x units per pack
Compatibility	ISO Class 4 and Grade A environments	
Compliance	BS EN ISO 23907, UN3291, AFNO NFX 30-511 and ISO 9002	
Sterilisation	Gamma irradiation	
Shelf Life	Sterile	3 years from the date of manufacture

\*triple bagged packaging available, enquire for a quotation

Helapet Limited | Lister House | Blackburn Road | Houghton Regis | Bedfordshire | LU5 5BQ

**Web:** [helapet.co.uk](http://helapet.co.uk) | **Phone:** +44 (0)1582 501980 | **Email:** [sales@helapet.co.uk](mailto:sales@helapet.co.uk)

Helapet Limited is a CliniMed Group Company.



**INVESTORS  
IN PEOPLE**  
Silver



CERTIFICATE NO. MD 78785