



berner

Cytotoxic Spill Kit

The Berner Cytotoxic Spill Kit contains a quality selection of compliant Personal Protective Equipment (PPE) and necessary aids to contain and remove hazardous cytotoxic spills.

Compliant with European PPE standards and CE marked, the Berner Cytotoxic Spill Kit enables quick and safe access by packing the contents in the order that the user will need them. All protective barrier items have proven permeation tests against a list of commonly used chemical and cytotoxic substances.

Each spill kit comes complete with full instructions for use, accident report form and handy illustrative guide.

Helapet supply the Berner Cytotoxic Spill Kit in packs of two, available in a range of four convenient garment configurations.



CS4004

Features and Benefits

- ✓ Contains essential PPE accessories
- ✓ Comprehensive permeation tests
- ✓ Proven barrier function
- ✓ Enclosed accident report form
- ✓ Gown and glove size options

Ordering Information

CS4004

Berner Cytotoxic Spill Kit with size 7.5 latex glove & medium gown non sterile • 2 pcs/pck

CS4004N

Berner Cytotoxic Spill Kit with size 7.5 neoprene glove & large gown non sterile • 2 pcs/pck

CS4004A

Berner Cytotoxic Spill Kit with size 9 latex glove & large gown non sterile • 2 pcs/pck

CS4004L


Berner Cytotoxic Spill Kit with size 9 neoprene glove & large gown non sterile • 2 pcs/pck



Specification

Usage	Single Use / Wet Spills Only	
Contents	1 x pair of cytotoxic protection gloves 1 x pair of covering gloves 1 x cytotoxic protection gown 1 x pair of shoe covers 1 x protective goggles 1 x particle-filtering P3 half-mask 2 x ChemoSorb pads 3 x cleaning cloths 1 x blue waste bag 1 x waste bag with imprint 2 x plastic bag locks 1 x instructions for use	
Glove sizes	Size 7.5 latex glove	CS4004
	Size 9 latex glove	CS4004A
	Size 7.5 neoprene glove	CS4004N
	Size 9 neoprene glove	CS4004L
Gown sizes	Medium Gown	CS4004
	Large Gown	CS4004A
	Large Gown	CS4004N
	Large Gown	CS4004L
Compliance	CE Marked; EC type tested and certified; ISO9001	
Shelf Life	3 years from date of manufacture	

Use of SpillKits for the disposal of spilled cytotoxins



Free 11 Point Procedure Poster

Request your free 11 point procedure poster on the proper emergency application of Berner Cytotoxic Spill Kits.

Enquire online or call our office today.

Helapet Limited

Lister House | Blackburn Road | Houghton Regis | Bedfordshire | LU5 5BQ

Web: helapet.co.uk | **Phone:** +44 (0)1582 501980 | **Email:** sales@helapet.co.uk

Helapet Limited is a CliniMed Group Company.

