



# berner

## Cytotoxic Workmat

Incorporating a proven 3 ply blue-mesh structure, the Berner Cytotoxic Workmat is designed for effective absorption and containment of accidental cytotoxic drips and spills.

A lint-free polypropylene top layer provides effective product protection during preparation activities, with an absorbent middle layer to contain and distribute spilt agents without losing shape. A non-slip bottom layer provides additional stability and prevents any spilt liquids from seeping onto the surface.

Available from Helapet in two sizes and suitable for use in latex-free areas.



### Features and Benefits

- ✓ Proven 3 layer protective system
- ✓ Extra absorbent
- ✓ Non-slip backing
- ✓ Reliable cytotoxic barrier properties
- ✓ Lint free and sterile


### Ordering Information

#### WP8001

26 x 42cm Berner Cytotoxic Workmat  
sterile • 50 pcs/pck

#### WP8101

41 x 56cm Berner Cytotoxic Workmat  
sterile • 40 pcs/pck



WP8101



## Technical Data

Product Code		WP8001	WP8101
Material	Under layer	Polyethylene foil	
	Central layer	Air-laid paper	
	Mesh layer	Polyethylene	
Dimensions W x L		42 x 26cm ±10%	56 x 41cm ±10%
Weight per box		16.5g	30g
Colour		Blue	
Absorption of liquid*		≤100ml	≤200ml
Packaging		50 pieces per dispenser carton	40 pieces per dispenser carton
Sterilisation		Ethylene Oxide (ETO)	
Compatibility		ISO Class 4 and Grade A environments	
Shelf life		4 years from date of manufacture	

\*per prep mat