



berner XP

Cytotoxic Spill Kit

The Berner XP Cytotoxic Spill Kit provides a complete emergency cleaning solution for wet and dry chemical and cytotoxic spills. Each Berner XP Spill Kit incorporates a quality list of latex-free clothing, tools and materials to safely manage and remove spilt medications and agents.

Compliant with European PPE standards and CE marked, the Berner XP Cytotoxic Spill Kit enables quick and safe access by packing the contents in the order that the user will need them. All protective barrier items have proven permeation tests against a list of commonly used chemical and cytotoxic substances.

Each spill kit comes complete with full instructions for use, accident report form and handy illustrative guide.



CS41XX

Features and Benefits

- ✓ Suitable for wet and dry spills
- ✓ Latex free
- ✓ Complete range of essential PPE accessories
- ✓ Proven barrier function
- ✓ Enclosed accident report form
- ✓ User instructions and illustration guide

Ordering Information

CS4124

Berner XP Cytotoxic Spill Kit
medium • non sterile • 1 pc/pck

CS4133

Berner XP Cytotoxic Spill Kit
large • non sterile • 1 pc/pck

CS4135

Berner XP Cytotoxic Spill Kit
XX large • non sterile • 1 pc/pck



Specification

Usage	Single Use / Wet and Dry Spills	
Contents	1 x protective Active Breath coverall 1 x pair of neoprene protective gloves 1 x pair of vinyl over gloves 1 x pair protective Active Breath overboots 1 x FFP3 respirator mask (latex-free elastomer bands) 1 x pair of safety goggles (latex-free elastomer bands) 1 x wax crayon 1 x warning sign 1 x ChemoSorb pad 4 x ISYSOFT cloths 1 x scoop 1 x sweeper 1 x tongs 1 x 500ml distilled water 1 x waste bag, blue 1 x special waste bag, white 2 x cable ties 1 x user information	
Glove sizes	Size 7.5 neoprene glove	CS4124
	Size 8.5 neoprene glove	CS4133
	Size 9 neoprene glove	CS4135
Gown sizes	Medium Gown	CS4124
	Large Gown	CS4133
	XX Large Gown	CS4135
Compliance	CE Marked; EC type tested and certified; ISO9001	
Shelf Life	3 years from date of manufacture	



Free 11 Point Procedure Poster

Request your free 11 point procedure poster on the proper emergency application of Berner Cytotoxic Spill Kits.
Enquire online or call our office today.

Helapet Limited

Lister House | Blackburn Road | Houghton Regis | Bedfordshire | LU5 5BQ

Web: helapet.co.uk | Phone: +44 (0)1582 501980 | Email: sales@helapet.co.uk

Helapet Limited is a CliniMed Group Company.

