



## Advance Cleanroom Glove

BioClean™ Advance Cleanroom Gloves offer exceptional comfort and dexterity, anatomically shaped for enhanced fit and easier grip. The robust latex material is compliant to EN ISO 374 standards, whilst maintaining exceptional flexibility. The smooth glove surface greatly assists double gloving where required.

BioClean™ Advance Cleanroom Gloves are cleanroom processed and supplied in paper-free packaging, ensuring compatibility up to ISO Class 4 & EU GMP Grade A cleanroom environments. A comprehensive range of gloves sizes are available to choose from.

### Features and Benefits

- ✓ Anatomically shaped and designed to be double-donnable
- ✓ Compatible up to ISO Class 4 and EU GMP Grade A
- ✓ Compliant to EN ISO 374
- ✓ Paper-free packaging
- ✓ Sterile

### Ordering Information

#### GL5065

BioClean™ Advance Cleanroom Glove 6.5 • sterile • 50 pairs/pck

#### GL5070

BioClean™ Advance Cleanroom Glove 7 • sterile • 50 pairs/pck

#### GL5075

BioClean™ Advance Cleanroom Glove 7.5 • sterile • 50 pairs/pck

#### GL5080

BioClean™ Advance Cleanroom Glove 8 • sterile • 50 pairs/pck

#### GL5085

BioClean™ Advance Cleanroom Glove 8.5 • sterile • 50 pairs/pck



Enquire today for  
**FREE**  
samples



### Technical Data

Material	Latex
Colour	Natural
Cuff Length	300mm
Packaging	Non-particulating EasyTear PE packaging 1 x pair per inner PE wallet; 1 x wallet per sealed PE pouch; 10 x pouches per sealed outer bag; 5 x outer bags per carton
Processing	Processed in an NEBB certified ISO Class 4 environment
AQL	0.65
Compatibility	ISO Class 4 & EU GMP Grade A
Compliance	Category 3 - Complex Design Personal Protective Equipment PPE Regulation (EU)2016/425 EN420:2003 +A1:2009 and EN374-1:2016, EN374-5:2016
Sterilisation	Gamma irradiation
Shelf Life	5 years from date of manufacture

BioClean™ is a registered trademark owned by Ansell Limited or one of its affiliates.

#### Helapet Limited

Lister House | Blackburn Road | Houghton Regis | Bedfordshire | LU5 5BQ

**Web:** [helapet.co.uk](http://helapet.co.uk) | **Phone:** +44 (0)1582 501980 | **Email:** [sales@helapet.co.uk](mailto:sales@helapet.co.uk)

Helapet Limited is a CliniMed Group Company.

