



Scan Me for Web Page



Safestep Non-slip Overshoes

Manufactured from lightweight, low linting polypropylene material, BioClean™ Safestep Non-slip Overshoes are highly durable and economic footwear with an elastomer coating to provide enhanced grip on smooth floors and surfaces.

Supplied in one standard size with an elasticated ankle, Safestep Non-slip Overshoes are designed to universally fit over different types of footwear. Each overshoe is low-particle shedding and cleanroom processed to ensure ISO Class 5 compatibility.

Features and Benefits

- ✓ Low linting materials
- ✓ Compatibility with ISO Class 5 conditions
- ✓ Anti-slip sole for added grip
- ✓ Elastic ankle cuff to fit different types of footwear

Ordering Information

AO0016

BioClean™ Safestep Non-slip Overshoes
406mm (16") • non sterile • 4 x 100 pcs/pck



Enquire today for **FREE** samples

AO0016



Technical Data

Material	Elastomer coating over non-woven spun-bonded polypropylene	
Colour	White	
Properties	Low particle shedding	
Material Weight	42 g/m ²	
Coating	Elastomer	
Thickness	0.14mm	
Dimensions	Height	170mm (6.5")
	Length	400mm (16")
Packaging	100 x pieces per sealed inner PE bag; 1 x inner bag per sealed outer PE bag; 4 x outer bags per lined carton	
Processed	NEBB certified ISO Class 5 cleanroom	
Compliance	ISO Class 5	
Sterilisation	Non sterile	
Shelf Life	5 years from date of manufacture	

BioClean™ is a registered trademark owned by Ansell Limited or one of its affiliates.

Helapet Limited | Lister House | Blackburn Road | Houghton Regis | Bedfordshire | LU5 5BQ

Web: helapet.co.uk | **Phone:** +44 (0)1582 501980 | **Email:** sales@helapet.co.uk

Helapet Limited is a CliniMed Group Company.

