



Scan Me for Web Page



# Combi-Cap Syringe Cap

Helapet Combi-Cap provides secure and versatile sterile closure of syringes, catheters and open ended IV devices using both a male and female luer lock connection. Non-toxic and non-pyrogenic, each cap is supplied individually blister packed and sterile.

Helapet Combi-Caps are microbiologically and structurally integrity tested with PharmaPack® sterile BD luer lock syringes, in accordance with the NHS Method Document: Protocols for the Integrity Testing of Devices\*.

## Features and Benefits

- ✓ Male and female luer lock closure
- ✓ Non-toxic and non-pyrogenic
- ✓ Allows alcohol disinfection with no leaching
- ✓ Sterile, individually blister packed
- ✓ CE marked



IV2102

## Ordering Information

### IV2102

Combi-Cap Syringe Cap  
sterile • 10 x 100 pcs/pck

## Technical Data

Material	Polypropylene (PP)
Weight	0.33g
Colour	Red
Packaging	Individually blister packed; 100 x blisters per box; 10 x boxes per pack
Sterilisation	Ethylene Oxide (ETO)
Compliance	Integrity tested w/ PharmaPack® sterile BD luer lock syringes in accordance with NHS Protocols for the Integrity Testing of Devices*; CE Marked
Shelf Life	5 years from date of manufacture

\*NHS Method: Protocols for the Integrity Testing of Devices; 2nd Edition, revised April 2013

Enquire today for **FREE** samples

Helapet Limited | Lister House | Blackburn Road | Houghton Regis | Bedfordshire | LU5 5BQ

**Web:** helapet.co.uk | **Phone:** +44 (0)1582 501980 | **Email:** sales@helapet.co.uk

Helapet Limited is a CliniMed Group Company. SteriClean® is a registered trademark of CliniMed (Holdings) Limited.