



Easycap®



Applied using a simple “push on” design, Easycap® provides tamper-evident capping of open access ports and luer valves.

Easycap® can be easily secured to the end of our Vial2Bag® IV Transfer System, preventing possible tampering or contamination of the medication port.

Helapet supply the Easycap® sterile blister packed for smooth aseptic transfer into cleanroom environments.

Features and Benefits

- ✓ Secure application
- ✓ Simple push on design
- ✓ ETO sterile and blister packaged
- ✓ CE marked



IV3100



Ordering Information

IV3100
Easycap®
sterile • 15 x 8 pcs/pck

Technical Data

Material	Polypropylene (PP)
Colour	Red
Packaging	8 x caps per individual sterile blister; 15 x blisters per pack
Sterilisation	Ethylene Oxide (ETO)
Compliance	CE marked
Shelf Life	3 years from date of manufacture

Enquire today for
FREE
samples