



Foam Beard Mask

Our fibre-free Foam Beard Mask combines effective particle/moisture control with comfort, not restricting speech or breathing, goggle/glasses steaming or moisture build-up on the face. Comfortably locking behind the ear using a loop fastening, each mask can protect the wearer from unwanted contact as well as restricting loose particle and moisture from entering the cleanroom.

Foam Beard Masks can snugly fit around the nose and mouth, making them ideal to be worn in conjunction with cleanroom suits with a built-in facemask or on their own. Excellent absorbency performance ensures retention of particles and moisture for long time periods.

Features and Benefits

- ✓ Effective particle / moisture control
- ✓ Comfortable, breathable fit across the mouth and nose
- ✓ Does not restrict breathing or speech
- ✓ Can be worn with other facial protection

Ordering Information

AF4402

Foam Beard Mask
std • non sterile • 10 pcs/pck

Enquire today for
FREE
samples

AF4402



Technical Data

Colour	White
Properties	Fibre-free
Particles	>.5µ Via LPC* <750/cm ²
TNVR's	DI H2O <.9g/m ²
Ions - count of extractable (PPM)	
Chloride	<5
Sodium	<15
Potassium	<2
Packaging	Individually poly-bagged; 10 x beard masks per pack
Compatibility	ISO Class 4 & EU GMP Grade A
Sterilisation	Non sterile
Shelf Life	5 years from date of manufacture

* Liquidborne Particle Counts

Helapet Limited

Lister House | Blackburn Road | Houghton Regis | Bedfordshire | LU5 5BQ

Web: helapet.co.uk | **Phone:** +44 (0)1582 501980 | **Email:** sales@helapet.co.uk

Helapet Limited is a CliniMed Group Company.

