



Glove Liner



The Helapet Glove Liner is a white, 18 gauge inner glove knitted with a highly stretchable mixture of nylon, acrylic fibres and stretchable polyamide. This material has an exceptional capacity to retain perspiration between its fibres, keeping the hand dry and free from moisture.

Completely lint-free, this glove liner can be ideally used as an underglove in cleanroom environments as well as a hand warmer for activities in colder working conditions.

The thin material grants easy double-gloving without losing touch sensitivity or dexterity.



Features and Benefits

- ✓ Lint free material
- ✓ Effective retention of perspiration and moisture
- ✓ Compatible as an underglove or hand warmer
- ✓ Grants easy double-gloving

Ordering Information

GW0900
 Glove Liner
 std • non sterile • 30 pairs/pck



Technical Data

Material	Vynylon - 70% Nylon, 30% Acrylic
Colour	White
Properties	100% Lint free; Retains perspiration; Washable at 40°C using a neutral detergent.
Packaging	10 x pairs per sealed polybag; 3 x polybags per pack
Compliance	PPE Category 1
Sterilisation	Non sterile
Shelf Life	5 years from date of manufacture