



# Intrathecal Bag

Transparent, disposable bag designed for secure site to site transportation of intrathecal drugs. The Helapet Intrathecal Bag is manufactured from strong, clear polyethylene material for easy visibility of contents, with a red and white colour coded front panel for easy labelling and logging of information.

Alcohol resistant ink allows liberal spraying with alcohol for disinfecting purposes. The sealable fold over closure allows for tamper-evident onetime use, with a perforated line for easy access. Available non sterile with full lot traceability.



## Features and Benefits

- ✓ Clear packaging to minimise errors
- ✓ Secure tamper-evident seal
- ✓ Perforated for easy access
- ✓ Tough 300 gauge / 75 micron polyethylene material
- ✓ Write-on panel for logging valuable information
- ✓ Alcohol resistant ink

PB0810

## Ordering Information

**PB0810**

W200 x L255mm • non sterile • 200 pcs/pck

Enquire today for  
**FREE**  
samples



## Technical Data

Dimensions (mm)	W200 x L255
Material	Polyethylene (PE)
Thickness	300 gauge / 75 micron
Packaging	200 x bags per resealable outer polybag
Sterilisation	Non sterile
Shelf Life	5 years from date of manufacture

Helapet Limited

Lister House | Blackburn Road | Houghton Regis | Bedfordshire | LU5 5BQ

**Web:** [helapet.co.uk](http://helapet.co.uk) | **Phone:** +44 (0)1582 501980 | **Email:** [sales@helapet.co.uk](mailto:sales@helapet.co.uk)

Helapet Limited is a CliniMed Group Company.

