



MedTrans Waste Bag

Helapet MedTrans Waste Bags are ideal for tamper-evident transportation of medications, or for the collection of waste in controlled working environments.

Manufactured from tough, clear polyethylene allows easy identification of the contents, with a write-on panel on one side for easy labelling. A tamper-evident fold over closure ensures secure containment of the product.

Sterile versions of the MedTrans Waste Bag are compatible for use in high grade cleanroom environments, including isolator units to package finished product or to collect waste.



PB1008

Features and Benefits

- ✓ Ideal for labelled transport or disposal of drugs
- ✓ Tough 400 gauge / 100 micron polyethylene material
- ✓ Secure tamper-evident seal
- ✓ Sterile and non sterile

Ordering Information

Sterile

XX1220

W320 x L450mm • sterile • 10 x 20 pcs/pck

Non sterile

PB1008

W320 x L450mm • non sterile • 4 x 50 pcs/pck

Enquire today for
FREE
samples



Technical Data

	XX1220	PB1220
Dimensions (mm)	W320 x L450	
Material	Polyethylene (PE)	
Thickness	400 gauge / 100 micron	
Packaging	20 x bags per easy-peel sterile pouch; 10 x pouches per sealed polybag; 1 x polybag per pack	50 x bags per resealable polybag; 4 x polybags per pack
Sterilisation	Gamma sterilised	Non sterile
Shelf Life	3 years from date of manufacture	5 years from date of manufacture

