



MiniPorter™

Carrier System

The ideal solution for transporting temperature sensitive samples, the Helapet MiniPorter™ is validated to securely accommodate a variety of temperature controlled products including vaccines and blood units.

Using fully replaceable components, the MiniPorter™ Carrier System can maintain the desired temperature threshold and provide effective thermal protection in variable weather conditions.

The robust outer bag incorporates a secure Velcro® fastening and long shoulder strap. Available in blue, a clear PVC window allows instructions and important documentation to be inserted and visible to users.



Features and Benefits

- ✓ Validated temperature control
- ✓ Includes Medicool® cool pack
- ✓ Secure Velcro® fastening
- ✓ Clear PVC Window for inserting delivery data
- ✓ Long over-the-shoulder strap

Ordering Information

Vaccine MiniPorter™ +2°C to +8°C **VPM011**

VPM011

Vaccine MiniPorter™

Includes 1 x Medicool® 11 and Medimat

260 x 179 x 131mm • non sterile • 1 pc/pck

VPI011

Refurbishment Pack Vaccine MiniPorter™

Includes replacement inner, 1 x Medicool® 11 and Medimat

Exclusive to Vaccine MiniPorter™ • non sterile • 1 pc/pck

Blood MiniPorter™ +2°C to +10°C **BPM011**

BPM011

Blood MiniPorter™

Includes 1 x Medicool® 11

260 x 179 x 131mm • non sterile • 1 pc/pck

BPI011

Refurbishment Pack Blood MiniPorter™

Includes replacement inner and 1 x Medicool® 11

Exclusive to Blood MiniPorter™ • non sterile • 1 pc/pck



VPM011



BPM011



Technical Data

Vaccine MiniPorter™

Colour	Blue
Temperature Threshold	+2°C to +8°C
Duration at Threshold Temperature	18 hours*
Product Volume (L)	0.73
Inner Moulding Composition	Polystyrene
Outer Case Composition	Nylon
Cold Source	1 x Medicoool® 11 and 1 x Medimat
External Dimensions (w x l x h)	260 x 179 x 131mm
Internal Dimensions (w x l x h)	210 x 115 x 30mm
Empty weight (kg)	0.73

Blood MiniPorter™

Colour	Blue
Temperature Threshold	+2°C to +10°C
Duration at Threshold Temperature	3.5 hours*
Product Volume (L)	0.73
Inner Moulding Composition	Polystyrene
Outer Case Composition	Nylon
Cold Source	1 x Medicoool® 11
External Dimensions (w x l x h)	260 x 179 x 131mm
Internal Dimensions (w x l x h)	210 x 130 x 75mm
Empty weight (kg)	0.73

* Validation information available on request

Helapet Limited

Lister House | Blackburn Road | Houghton Regis | Bedfordshire | LU5 5BQ

Web: helapet.co.uk | Phone: +44 (0)1582 501980 | Email: sales@helapet.co.uk

Helapet Limited is a CliniMed Group Company.

