




Scan Me for Web Page

# Polypropylene Containment Socks



Designed to absorb large quantities of aggressive liquid leaks and spills, Helapet Polypropylene Containment Socks are ideal for use in cleanrooms, laboratories or other controlled environments where potential chemical containment and absorption is required.

Polypropylene Containment Socks prevent the spreading of liquid substances by easily shaping around the spill, reducing the contaminated area and allowing for quick and efficient clean-up. Each sock is individually roll packed in a sealed polybag, ensuring the whole product is free from unwanted debris.

CS2004



## Features and Benefits

- ✓ Low-linting polypropylene materials
- ✓ Can circle spills preventing migration
- ✓ Suitable for soaking chemical spills
- ✓ Rolled in heat sealed polybags

Enquire today for **FREE** samples

## Ordering Information

### CS2004

Polypropylene Containment Socks non sterile • 12 pcs/pck

## Technical Data

Material	Skin	Polypropylene
	Filler	Polypropylene Pulp
Dimensions	W8 x L117cm	
Size of outer bag	25cm x 38cm	
Weight	5.08kg	
Colour	Pink	
Packaging	Individually poly bagged and heat sealed; 12 x socks per pack	
Compatibility	Acids, Caustics and unknowns	