



Polythene Oversleeves

Polythene Oversleeves offer cost effective splash resistant protection to arms and garment sleeves. The lightweight, low linting material provides enhanced arm protection with a universal fit over standard cleanroom suits and workwear. Elasticated edges secure the oversleeve in place without restricting freedom of movement.

Helapet offer our Polythene Oversleeves non sterile or terminally sterilised in pairs.

Features and Benefits

- ✓ Splash resistant protection
- ✓ Light and comfortable
- ✓ Low linting materials
- ✓ Sterile and non sterile
- ✓ 40cm length

Ordering Information

AN0068

Polythene Oversleeves
std • non sterile • 100 pairs/pck

XX0068

Polythene Oversleeves
std • sterile • 100 pairs/pck



Enquire today for
FREE
samples

AN0068/ XX0068



Technical Data

Material	Low-density polyethylene (LDPE)	
Colour	White	
Length	40cm	
Properties	Concealed elasticated wrist	
Thickness	18 - 20 microns	
Packaging	Sterile	1 x pair per sealed polybag, 10 x bags per outer re- sealable bag; 10 x re-sealable bags per outer polybag.
	Non sterile	50 x pairs per sealed polybag; 2 x polybags per carton
Sterilisation	Gamma irradiation	
Shelf Life	Sterile	3 years
	Non sterile	5 years

Helapet Limited

Lister House | Blackburn Road | Houghton Regis | Bedfordshire | LU5 5BQ

Web: helapet.co.uk | **Phone:** +44 (0)1582 501980 | **Email:** sales@helapet.co.uk

Helapet Limited is a CliniMed Group Company.

