

Prep Trays

Sterile and Non sterile



Helapet offer a quality range of sterile and non-sterile Prep Trays. Durable and highly versatile, Helapet Prep Trays support safe handling of medical devices, sharps and medications, facilitating transfer of materials in and out of isolator units.

Available in two convenient sizes, sterile Prep Trays are individually polybagged* and gamma sterilised at no less than 25 kGy.

*Double and triple-bagged packaging available on request. Enquire for a quotation.

Features and Benefits

- ✓ Lightweight, durable design
- ✓ Colour: blue
- ✓ Easy to clean
- ✓ Two size options
- ✓ Sterile and non sterile



PT1400

Ordering Information



Sterile

XX2312

23 x 12 x 2cm • sterile • 2 x 50 pcs/pck

XX2428

28 x 23 x 5 cm • sterile • 2 x 50 pcs/pck

Non sterile

PP2312

23 x 12 x 2cm • non sterile • 30 x 5 pcs/pck

PP2428

28 x 23 x 5 cm • non sterile • 15 x 5 pcs/pck





Technical Data

Sterile

	XX2312	XX2428
Dimensions (l x w x h)	230 x 120 x 20mm	280 x 230 x 50mm
Colour	Sky blue	
Packaging	Individually polybagged; 25 x polybags per sealed outer bag, 2 x outer bags per pack	
Compliance	Suitable for use in ISO Class 4 and Grade A cleanrooms	
Processing	Individually cleaned and packaged in an ISO Class 7 cleanroom	
Sterilisation	Gamma irradiated	
Shelf Life	3 years from the date of manufacture	

Non sterile

	PP2312	PP2428
Dimensions (l x w x h)	230 x 120 x 20mm	280 x 230 x 50mm
Colour	Sky blue	
Packaging	5 x trays per heat sealed polybag; 30 x polybags per pack	5 x trays per heat sealed polybag; 15 x polybags per pack
Compliance	Suitable for use in ISO Class 4 and Grade A cleanrooms	
Processing	Individually cleaned and packaged in an ISO Class 7 cleanroom	
Sterilisation	Non sterile	
Shelf Life	5 years from the date of manufacture	

