



# Sterile Resealable Polybag

Helapet offer a range of sterile ready-to-use resealable polybags, ideal for bagging and labelling of small sterile samples to larger products.

Available in three size options, the clear polyethylene material allows for easy identification of the contents with a resealable tab for multi-opening. A4 size options can comfortably hold worksheets and other important paperwork for convenient transfer into and between cleanroom environments.

Each bag is packaged in an ISO Class 4 environment, ready folded and placed in a heat-sealed easy peel pouch limiting the risk of dust or moisture contamination. Each bag is gamma sterilised to no less than 25 kGy.



## Features and Benefits

- ✓ Robust 250 gauge / 62 micron polyethylene material
- ✓ Comprehensive range of sizes including A4
- ✓ Resealable for multi-use activities
- ✓ Heat sealed and packed in an easy-peel pouch
- ✓ Sterile

## Ordering Information

### XX1004

W228 x L323mm (A4) • sterile • 8 x 25 pcs/pack

### XX1005

W228 x L323mm (A4) • sterile • 20 x 10 pcs/pack

### XX1116

W279 x L406mm • sterile • 40 x 5 pcs/pack

### XX1117

W279 x L406mm • sterile • 20 x 10 pcs/pack

### XX1375

W76 x L190mm • sterile • 40 x 15 pcs/pck

Enquire today for  
**FREE**  
samples



## Technical Data

	<b>XX1004</b>	<b>XX1005</b>	<b>XX1116</b>	<b>XX1117</b>	<b>XX1375</b>
Dimensions (mm)	W228 x L323 (A4 paper)	W228 x L323 (A4 paper)	W279 x L406	W279 x L406	W76 x L190
Material	Virgin Low Density Polyethylene (LDPE)				
Thickness	250 gauge / 62 micron				
Packaging	5 x bags per easy-peel sterile pouch; 5 x pouches per re-sealable polybag; 8 x bags per pack	10 x bags per easy-peel sterile pouch; 20 x pouches per sealed polybag; 1 x polybag per pack	5 x bags per easy-peel sterile pouch; 40 x pouches per sealed polybag; 1 x polybag per pack	10 x bags per easy-peel sterile pouch; 20 x pouches per sealed polybag; 1 x polybag per pack	15 x bags per easy-peel sterile pouch; 40 x pouches per sealed polybag; 1 x polybag per pack
Sterilisation	Gamma sterilised				
Shelf Life	3 years from date of manufacture				

Helapet Limited

Lister House | Blackburn Road | Houghton Regis | Bedfordshire | LU5 5BQ

**Web:** helapet.co.uk | **Phone:** +44 (0)1582 501980 | **Email:** sales@helapet.co.uk

Helapet Limited is a CliniMed Group Company.

