

Scan Me for Web Page



Sempermed[®] Supreme Surgical Glove

Helapet offer a comprehensive range of SemperMed® Supreme latex gloves. Each glove is anatomically shaped for precise handling and safety during surgical procedures.

A fully synthetic inner coating allows for simple and easy donning, with a micro-rough outer surface for secure grip. The thinness of the material allows for double-gloving without reducing sensitivity. A long 300mm beaded cuff provides a secure fit and additional protection above the wrist.

SemperMed® Supreme Surgical Gloves are cleanroom processed in paper-free packaging and presented gamma sterile. Compatible up to ISO Class 4 & EU GMP Grade A operating conditions.

Features and Benefits

- Anatomically shaped
- Easy donning and double gloving
- ✓ Compatible up to ISO Class 4 & EU GMP Grade A
- ✓ Micro-rough surface for enhanced grip
- ✓ Compliant to EN374
- ✓ Sterile

Ordering Information

GL1060-50

Sempermed® Supreme Surgical Glove 6 • sterile • 50 pairs/pck

GI 1065-50

Sempermed® Supreme Surgical Glove 6.5 • sterile • 50 pairs/pck

GL1070-50

Sempermed® Supreme Surgical Glove 7 • sterile • 50 pairs/pck

GL1075-50

Sempermed® Supreme Surgical Glove 7.5 • sterile • 50 pairs/pck

GL1080-50

Sempermed® Supreme Surgical Glove 8 • sterile • 50 pairs/pck

GL1085-50

Sempermed® Supreme Surgical Glove 8.5 • sterile • 50 pairs/pck

GL1090-50

Sempermed® Supreme Surgical Glove 9 • sterile • 50 pairs/pck



Helapet Limited | Lister House | Blackburn Road | Houghton Regis | Bedfordshire | LU5 5BQ Web: helapet.co.uk | Phone: +44 (0)1582 501980 | Email: sales@helapet.co.uk







Technical Data

Material	Natural Rubber Latex
Colour	White
Cuff Length	300mm
Properties	Anatomical with beaded cuff
Packaging	1 x pair in sealed medical peel pack; 50 x peel packs per dispenser
Processing	Processed in an ISO Class 4 environment
AQL	0.65
Compatibility	ISO Class 4 & EU GMP Grade A operating conditions
Compliance	Category 3 PPE: 89/686/EEC
	EN-374-2, EN 1041, EN 980, EN 455/3, ISO 15223
Sterilisation	Gamma irradiation
Shelf Life	3 years from manufacture

