



Steri-Dropper™ Bottles

Steri-Dropper™ Sterile Dropper Bottles are a latex free, non-toxic and non-pyrogenic alternative to the traditional glass dropper bottle.

Avoiding the need for lengthy and expensive sterilisation measures, Steri-Dropper™ Bottles come in sterile packs and are ideal for accurate preparation and dispensing of sterile formulations in cleanroom environments.

Provided in a range of four sizes, each pack contains two separate resin free bottles with snap-in dropper tips and polypropylene screw caps.

Supplied sterile in an easy peel blister with low-linting Tyvek® backing.

Features and Benefits

- ✓ Easy squeeze polyethylene
- ✓ Natural colour for easy visual identification
- ✓ Available in four different sizes
- ✓ Recyclable blister packaging
- ✓ Sterile packaged in pairs



Ordering Information

IV1630

3ml • sterile • 50 x 2 pcs/pck

IV1634

10ml • sterile • 50 x 2 pcs/pck

IV1632

7ml • sterile • 50 x 2 pcs/pck

IV1633

15ml • sterile • 50 x 2 pcs/pck

Enquire
today for
FREE
samples

Technical Data

Steri-Dropper™ Bottles

	IV1630	IV1632	IV1634	IV1633
Bottle size	3ml	7ml	10ml	15ml
Colour	Clear, natural			
Material	Low density polyethylene and zinc stearate-free resin, latex free			
Packaging material	Clear plastic moulded blister, with Tyvek peel away lid			
Packaging	2 x bottles per easy peel sterile blister w/ Tyvek® lid; 50 x sterile blisters per pack			
Sterilisation	Ethylene Oxide (ETO)			
Shelf life	2 years			

Helapet Limited

Lister House | Blackburn Road | Houghton Regis | Bedfordshire | LU5 5BQ

Web: helapet.co.uk | **Phone:** +44 (0)1582 501980 | **Email:** sales@helapet.co.uk

Helapet Limited is a CliniMed Group Company.

