

# Steri-Dropper™ Sterile Eye-Dropper Bottles



Steri-Dropper™ Steri-Dropper Sterile Eye-Dropper Bottles are latex free and a non-toxic, non-pyrogenic alternative to the traditional glass eye-dropper bottle.

Avoiding the need for lengthy and expensive sterilisation measures, Steri-Dropper™ Bottles come in sterile packs and are ideal for accurate dispensing of formulations in cleanroom environments.

Provided in a range of four sizes, each pack contains two separate resin free bottles with snap-in dropper tips and polypropylene screw caps. Supplied sterile in an easy peel blister with low-linting Tyvek® backing.

## Features and Benefits

- ✓ Easy squeeze polyethylene
- ✓ Natural colour for easy visual identification
- ✓ Available in four different sizes
- ✓ Recyclable blister packaging
- ✓ Sterile



## Ordering Information

### IV1630

3ml • sterile • 50 x 2 pcs/pck

### IV1634

10ml • sterile • 50 x 2 pcs/pck

### IV1632

7ml • sterile • 50 x 2 pcs/pck

### IV1633

15ml • sterile • 50 x 2 pcs/pck

Enquire today for  
**FREE**  
samples

## Technical Data

**Steri-Dropper™ Sterile Eye Dropper Bottles**

	<b>IV1630</b>	<b>IV1632</b>	<b>IV1634</b>	<b>IV1633</b>
Bottle size	3ml	7ml	10ml	15ml
Colour	Clear, natural			
Material	Low density polyethylene and zinc stearate-free resin, latex free			
Packaging material	Clear plastic moulded blister, with Tyvek peel away lid			
Packaging	2 x bottles per easy peel sterile blister w/ Tyvek lid; 50 x sterile blisters per pack			
Sterilisation	Ethylene Oxide (ETO)			
Shelf life	2 years			

