



SteriClean®

Alcohol Sprays

High performance general purpose disinfectant sprays for use in cleanrooms and critically controlled environments. SteriClean® Alcohol Sprays use an optimal mix of 70% v/v isopropyl alcohol (IPA) or de-natured ethanol (DE), blended with 30% water for injection (WFI) and 0.2 micron filtered for use in aseptic areas.

Contained in an ergonomic and lightweight bottle, each spray uses a proven closed bag in bottle spray assembly to maintain sterility of contents by restricting the ingress of bacteria*. 0.45L trigger sprays are ideal for accessible cleaning in confined work areas. Each SteriClean® Alcohol Spray is double or triple bagged in an ISO Class 7 cleanroom and gamma sterilised to no less than 25 kGy, suitable for efficient aseptic transfer into cleanrooms and isolator environments.

Features and Benefits

- ✓ 70/30 isopropyl alcohol (IPA) or denatured ethanol (DE)
- ✓ 0.2 micron filtered
- ✓ Bag in bottle closed system ensures sterility of contents
- ✓ Lightweight, ergonomic bottle design
- ✓ Long length, non-venting trigger
- ✓ 0.9L and 0.45L[†] bottle sizes
- ✓ Individually double or triple bagged and gamma sterilised

*Data available on request
†0.45L IPA Spray - minimum order quantity of 180 pcks



Triple bagged packaging available
Enquire for a quotation

Enquire today for
FREE
samples



XX6800

XX6902



Ordering Information

0.9L Bottles

XX6800

SteriClean® IPA Spray (bag in bottle)
0.9L • sterile • 6 bottles/pck

XX6900

SteriClean® DE Spray (bag in bottle)
0.9L • sterile • 6 bottles/pck

XXTB6900

SteriClean® DE Spray (bag in bottle)
triple bagged
0.9L • sterile • 6 bottles/pck

XX6901

SteriClean® DE Spray (no bag in bottle)
0.9L • sterile • 6 bottles/pck

SP6801

SteriClean® IPA Spray (no bag in bottle)
0.9L • non sterile • 6 bottles/pck

0.45L Bottles

XX6902

SteriClean® DE Spray (bag in bottle)
0.45L • sterile • 12 bottles/pck



SP6801

Technical Data

Solution	70% Isopropyl Alcohol (IPA) or Denatured Ethanol (DE)	
	30% Water for Injection	
	0.2 µm filtered	
Bottle	Colour	White
	Properties	Bag in bottle closed system and non-venting trigger
Size	0.9L and 0.45L	
Processing	Manufactured and packaged in ISO Class 7 cleanroom	
	Filled in ISO Class 5 environment	
Packaging	Individually double or triple polybagged	
Sterility	Sterile products are gamma irradiated to no less than 25kGy	
Validated Shelf Life	2 years from date of manufacture	
In-use Shelf Life	3 months from date of opening*	

*Data available on request

Microbiological Efficacy

Test	Substance	Minimum Contact Time
EN 1276	Bacteria	5mins
EN 13697	Bacteria	5mins
EN 13697	Fungi	15mins
EN 1650	Fungi	15mins

Helapet Limited

Lister House | Blackburn Road | Houghton Regis | Bedfordshire | LU5 5BQ

Web: helapet.co.uk | Phone: +44 (0)1582 501980 | Email: sales@helapet.co.uk

Helapet Limited is a CliniMed Group Company.

