

SteriClean®

Dry Wipes



SteriClean® Dry Wipes are made from polyester/cellulose bonded fibres which form an extremely robust and absorbent material.

Low linting with low extractables, SteriClean® Dry Wipes are supplied ready folded in alcohol resistant paper free packaging and irradiated at not less than 25 kGy making them an ideal choice for ISO class 5 cleanrooms.

Helapet supply a larger 30cm x 30cm wipe, ideal for economical cleaning of larger surface areas with fewer wipes.

Available in several packaging configurations from individually packaged to economic multi-packs, all batches are cleanroom processed with lot traceability.

Features and Benefits

- ✓ 68gsm polyester / cellulose material
- ✓ Absorbent, low linting with low extractables
- ✓ Two wipe size options
- ✓ Cleanroom processed
- ✓ Packed in alcohol resistant foil pouch
- ✓ Sterile and non sterile presentations



Enquire today for **FREE** samples

Ordering Information

Standard Dry Wipes

XX6700

SteriClean® Single Pack Dry Wipe
23 x 23cm • sterile • 300 x 1pc/pck

XX6705

SteriClean® 5 Pack Dry Wipe
23 x 23cm • sterile • 60 x 5 pcs/pck

XX6710

SteriClean® 10 Pack Dry Wipe
23 x 23cm • sterile • 30 x 10 pcs/pck

XX6709

SteriClean® (Bulk Pack) Dry Wipes
23 x 23cm • sterile • 3 x 300 pcs/pck

XL Dry Wipes

XX6712

SteriClean® XL 5 Pack Dry Wipes
30 x 30cm • sterile • 60 x 5 pcs/pck

XX6769

SteriClean® XL Dry Wipes
30 x 30cm • sterile • 10 x 50 pcs/pck

Non Sterile Dry Wipes

WN6709

SteriClean® Dry Wipes
23 x 23cm • non sterile • 300 pcs/pck

WN6712

SteriClean® XL Dry Wipes
30 x 30cm • non sterile • 150 pcs/pck



Technical Data

Wipe Size (unfolded)	Standard Dry Wipe	23 x 23cm
	XL Dry Wipe	30 x 30cm
Wipe Material	Low-linting 55% polyester / 45% cellulose Binder-free	
Colour	Natural	
Thickness	68gsm	
Absorbency	500 g/m ²	
Packaging	Standard and XL Dry Wipes	Alcohol resistant foil pouch with flat top seal
	Bulk Packed Dry Wipes	Double bagged; 2 x sealed 200g polyethylene bags
Processing	Ready folded and packaged in ISO Class 7 cleanroom	
Sterilisation	Sterile products are gamma irradiated to no less than 25kGy	
Shelf Life	Sterile	2 years from date of manufacture
	Non sterile	5 years from date of manufacture

Helapet Limited

Lister House | Blackburn Road | Houghton Regis | Bedfordshire | LU5 5BQ

Web: helapet.co.uk | **Phone:** +44 (0)1582 501980 | **Email:** sales@helapet.co.uk

Helapet Limited is a CliniMed Group Company.

