

SteriClean®

Foam Wipe

Produced from thick, soft polyurethane foam, SteriClean® Foam Wipes incorporate sponge-like characteristics which give excellent absorbency rates while still being suitable for use cleanrooms. Effective against accidental spills and soaking up large volumes of liquid, each wipe uses an open-cell foam structure for use on scratch sensitive surfaces and electronics. Very low linting and presented sterile, SteriClean® Foam Wipes are ideal for extra absorbency requirements and robust cleaning activities in aseptic working environments.


Features and Benefits

- ✓ High absorbency performance
- ✓ Excellent abrasion resistance
- ✓ Soft and easy to handle
- ✓ Low linting
- ✓ Ideal for spill control and cleaning of medical devices
- ✓ Individually poly-bagged
- ✓ Double bagged and sterile

Ordering Information

XX4699

SteriClean® Foam Wipe
23 x 23cm • sterile • 10 x 10 pcs/pck



XX4699

Enquire today for
FREE
samples



Technical Data

Wipe Size (unfolded)	23 x 23cm
Wipe material	Open cell, fibre-free foam construction
Colour	White
Thickness	0.109"
Particle >0.5 µm / m ²	<10,000,000/m ²
Packaging	Individually poly-bagged; 10 x polybags in outer bag; 10 x outer bags per pack
Processing	Packaged in ISO Class 7 cleanroom
Sterilisation	Sterile products are gamma irradiated to no less than 25kGy
Shelf Life	3 years from date of manufacture

