

SteriClean®

Individual IPA Wipes



Single sterile sachet wipes impregnated with our high performance 70% v/v IPA solution. SteriClean® Individual Wipes are perfect for precise cleaning and disinfection of medical products, devices and surfaces within cleanroom and isolator environments. Individually packed in paper-free foil packaging helps to optimise space in confined working environments and reduces waste compared to larger multi-packs. Each wipe is supplied pre-folded and sterile, available in both a standard and larger sizes.

Features and Benefits

- ✓ Single low linting, 68 gsm wipes
- ✓ High performance polyester / cellulose material
- ✓ Excellent economy cleaning and disposal
- ✓ Cleanroom processed and sterile
- ✓ Nick for easy opening
- ✓ Available in two convenient sizes



Ordering Information

XX6720

Individual IPA Wipe - Standard
10 x 15cm • sterile • 5 x 100 pcs/pck

XXTB6725

Individual IPA Wipe - Large
triple bagged
19 x 20cm • sterile • 10 x 25 pcs/pck



25 x Foil Sachets
+ 2 x Sealed Polybags

Enquire today for
FREE
samples



Technical Data

Solution	70% Isopropyl Alcohol (IPA)	
	30% Water for Injection (WFI)	
	0.2 µm filtered	
Wipe Size (unfolded)	XX6720 (Std)	10 x 15cm
	XXTB6725 (Lge)	19 x 20cm
Sachet Size	XX6720 (Std)	6 x 8cm
	XXTB6725 (Lge)	6 x 13cm
Wipe material	Low-linting polyester/cellulose	
Thickness	68 gsm	
Packaging	XX6720 (Std)	100 x foil sachets + 1 x resealable polybag; 5 x polybags per pack
	XXTB6725 (Lge)	25 x foil sachets + 2 x sealed polybags; 10 polybags per pack
Sterilisation	Sterile products are gamma irradiated to no less than 25kGy	
Validated Shelf Life	2 years from date of manufacture	

Microbiological Efficacy

Test	Substance	Minimum Contact Time
EN 1276	Bacteria	5mins
EN 13697	Bacteria	5mins
EN 13697	Fungi	15mins
EN 1650	Fungi	15mins

