

SteriClean®

Individual IPA Wipes



Single sterile sachet wipes impregnated with our high performance 70% v/v IPA solution. SteriClean® Individual Wipes are perfect for precise cleaning and disinfection of medical products, devices and surfaces within cleanroom and isolator environments. Individually packed in paper-free foil packaging helps to optimise space in confined working environments and reduces waste compared to larger multi-packs. Each wipe is supplied pre-folded and sterile, available in both a standard and larger sizes.

Features and Benefits

- ✓ Single low linting, 68 gsm wipes
- ✓ High performance polyester / cellulose material
- ✓ Excellent economy cleaning and disposal
- ✓ Cleanroom processed and sterile
- ✓ Nick for easy opening
- ✓ Available in two convenient sizes



Ordering Information

XX6720

Individual IPA Wipe - Standard
10 x 15cm • sterile • 5 x 100 pcs/pck

XXTB6725

Individual IPA Wipe - Large
triple bagged
19 x 20cm • sterile • 10 x 25 pcs/pck



25 x Foil Sachets
+ 2 x Sealed Polybags





Technical Data

Solution	70% Isopropyl Alcohol (IPA)	
	30% Water for Injection (WFI)	
	0.2 µm filtered	
Wipe Size (unfolded)	XX6720 (Std)	10 x 15cm
	XXTB6725 (Lge)	19 x 20cm
Sachet Size	XX6720 (Std)	6 x 8cm
	XXTB6725 (Lge)	6 x 13cm
Wipe material	Low-linting polyester/cellulose	
Thickness	68 gsm	
Packaging	XX6720 (Std)	100 x foil sachets + 1 x resealable polybag; 5 x polybags per pack
	XXTB6725 (Lge)	25 x foil sachets + 2 x sealed polybags; 10 polybags per pack
Sterilisation	Sterile products are gamma irradiated to no less than 25kGy	
Validated Shelf Life	3 years from date of manufacture	

Microbiological Efficacy

Test	Substance	Minimum Contact Time
EN 1276	Bacteria	5mins
EN 13697	Bacteria	5mins
EN 13697	Fungi	15mins
EN 1650	Fungi	15mins

Helapet Limited

Lister House | Blackburn Road | Houghton Regis | Bedfordshire | LU5 5BQ

Web: helapet.co.uk | **Phone:** +44 (0)1582 501980 | **Email:** sales@helapet.co.uk

Helapet Limited is a CliniMed Group Company.

