

SteriClean[®]

Laser Sealed Polywipe



Helapet offer a durable low-linting cleanroom wipe with excellent chemical compatibility. Laser sealed edges help to reduce the risk of loose fibres, essential for cleanroom environments.

Each polywipe is constructed from highly absorbent textured polyester knitted fabric, effective at lifting particles without leaving abrasions.

Features and Benefits

- ✓ Low linting laser-sealed edges for controlled environment cleaning
- ✓ Thin, durable PES material
- ✓ High absorbency of particles
- ✓ Binder-free
- ✓ Non sterile presentation

Ordering Information

WL6259
Laser Sealed Polywipe
23 x 23cm • non sterile • 150 pcs/pck

Enquire today for **FREE** samples

WL6259

Technical Data

Wipe Size (unfolded)	23 x 23cm
Wipe Material	100% interlock knitted polyester Binder free
Colour	White
Thickness	115gsm
Packaging	150 wipes per 1 x sealed polybag; 1 x polybag per pack
Processing	Packaged in an ISO Class 7 cleanroom environment