



# Tyvek®

## All-in-one Suit

Helapet offer our single-use Tyvek® All-in-one Suit for total body, head and foot coverage in high and low level cleanrooms.

Each suit offers superior levels of protection against liquids, chemical splashes and dry particles. Low-linting and abrasion resistant makes this suit ideal for robust use in controlled environments, whilst offering excellent breathability and comfort over long periods.

Available in a comprehensive range of sizes, the Tyvek® All-in-one Suit comes with an elasticated hood, cuffs and waist for optimum fit. A zip front closure with covering flap provides additional protection against snagging and possible contamination. Helapet offer this product non-sterile or double-bagged and terminally sterilised for higher grade cleanroom requirements.



### Features and Benefits

- ✓ Category III protective clothing
- ✓ Includes hood and overshoe protection
- ✓ Zip front closure with covering flap
- ✓ Low particle shedding
- ✓ Individually double-bagged and sterile, or non sterile presentation

### Ordering Information

#### Non sterile

##### AT0110

Tyvek® All-in-one Suit  
Small • non sterile • 25 pcs/pck

##### AT0120

Tyvek® All-in-one Suit  
Medium • non sterile • 25 pcs/pck

##### AT0130

Tyvek® All-in-one Suit  
Large • non sterile • 25 pcs/pck

##### AT0140

Tyvek® All-in-one Suit  
X Large • non sterile • 25 pcs/pck

##### AT0150

Tyvek® All-in-one Suit  
XX Large • non sterile • 25 pcs/pck

#### Sterile

##### XX0110

Tyvek® All-in-one Suit  
Small • sterile • 10 pcs/pck

##### XX0120

Tyvek® All-in-one Suit  
Medium • sterile • 10 pcs/pck

##### XX0130

Tyvek® All-in-one Suit  
Large • sterile • 10 pcs/pck

##### XX0140

Tyvek® All-in-one Suit  
X Large • sterile • 10 pcs/pck

##### XX0150

Tyvek® All-in-one Suit  
XX Large • sterile • 10 pcs/pck



Enquire today for  
**FREE**  
samples




## Technical Data

Material	Tyvek® flashspun polyethylene	
Colour	White	
Properties	Latex-free; low-linting; low particle generation	
Material Weight	180g per garment	
Packaging	Sterile	Individually double-bagged; 10 x sterile bags per pack
	Non sterile	5 x gowns per polybag; 5 x polybags per carton
Compliance	Category 3 - Complex Design Personal Protective Equipment PPE Directive 89/686/EEC	
Sterilisation	Ethylene Oxide	
Shelf Life	Sterile	3 years from date of manufacture
	Non sterile	5 years from date of manufacture

## Size Chart

Dimensions in centimetres

Size	A Chest	C Neck to Wrist	D Neck to Ankle	E Inside Leg
S	110	83	163	80
M	120	85	169	82
L	130	89	175	84
XL	140	92	183	86
XXL	150	96	192	88



Helapet Limited

Lister House | Blackburn Road | Houghton Regis | Bedfordshire | LU5 5BQ

Web: [helapet.co.uk](http://helapet.co.uk) | Phone: +44 (0)1582 501980 | Email: [sales@helapet.co.uk](mailto:sales@helapet.co.uk)

Helapet Limited is a CliniMed Group Company.

