



Tyvek®

Lab Coat

Helapet Tyvek® Lab Coats are created from lightweight non-woven spun bonded polyethylene material. Ideal for single-use partial bodily protection in less critical environments, Tyvek® Lab Coats can repel liquid splashes and dry particles.

Ideal for guests and staff, each lab coat offers elasticated cuffs for secure fit and protection around the wrist, fastened together with a frontal covered zip.

Helapet offer our Tyvek® Lab Coat in a range of sizes, with or without pockets.



Features and Benefits

- ✓ Low linting polyethylene Tyvek® material
- ✓ Available with or without pockets
- ✓ Protection from liquid splashes and dry particles
- ✓ Ideal for visitors

Ordering Information

No Pockets

AT1410

Tyvek® Lab Coat no pockets
Small • non sterile • 50 pcs/pck

AT1420

Tyvek® Lab Coat no pockets
Medium • non sterile • 50 pcs/pck

AT1430

Tyvek® Lab Coat no pockets
Large • non sterile • 50 pcs/pck

Two Pockets

AT1520

Tyvek® Lab Coat two pockets
Medium • non sterile • 50 pcs/pck

AT1530

Tyvek® Lab Coat two pockets
Large • non sterile • 50 pcs/pck

Enquire today for
FREE
samples





Technical Data

Material	Tyvek® flashspun polyethylene
Colour	White
Properties	Stitched internal seams, zip closure, elasticated cuffs, ultra-low linting, anti-statically treated
Packaging	25 x coats per outer bag; 2 x bags per carton
Type	No pockets; two pockets
Compliance	Category I PPE
	Antistatic treatment (EN 1149-5) - on both sides
Sterilisation	Non sterile
Shelf Life	5 years from date of manufacture

Size Chart

Dimensions in centimetres

Size	A Chest	B Back Length	C Neck to Sleeve	D Zip
S	108	107	77	70
M	118	109	78	70
L	128	111	79	70



Helapet Limited

Lister House | Blackburn Road | Houghton Regis | Bedfordshire | LU5 5BQ

Web: helapet.co.uk | Phone: +44 (0)1582 501980 | Email: sales@helapet.co.uk

Helapet Limited is a CliniMed Group Company.

