



Tyvek® Oversleeves



Tyvek® Oversleeves are composed of non-woven high density polyethylene to ensure the right combination of durability, resistance and comfort.

A proven non-linting and abrasion resistant material Tyvek® Oversleeves offer enhanced barrier protection against liquids, aerosols and fine particles when used in combination with chemical protective clothing.

The universal shape grants an easy and secure fit, able to be placed over other cleanroom apparel. Elasticated cuffs on both ends allow for a secure fit around the wrist and arm without restricting freedom of movement. Helapet supply Tyvek® Oversleeves in pairs, non-sterile or sterile wrapped in a low lint, sealed polybag.



Features and Benefits

- ✓ Durable non-woven fabric
- ✓ 42cm length
- ✓ Elasticated cuffs
- ✓ Sterile and non sterile

Ordering Information

AT0006

Tyvek® Oversleeves
non sterile • 100 pairs/pck

XX0006

Tyvek® Oversleeves
sterile • 100 pairs/pck

Enquire
today for
FREE
samples



Technical Data

| | | |
|---------------|--|---|
| Material | Tyvek® flashspun polyethylene | |
| Colour | White | |
| Length | 42cm | |
| Properties | Stitched internal seam; elasticated openings | |
| Packaging | Sterile | 1 x pair in a sealed polybag; 10 x bags per outer re-sealable bag; 10 x re-sealable bags per outer polybag. |
| | Non sterile | 50 x pairs per polybag, 2 x polybags per pack |
| Compliance | PPE Category I; Antistatic treatment (EN 1149-5) | |
| Sterilisation | Gamma irradiation | |
| Shelf Life | Sterile | 3 years from date of manufacture |
| | Non sterile | 5 years from date of manufacture |

