



UV Light Inhibiting Amber Bags

Certain medications, antibiotics and compounded chemotherapy products need to be stored and transported in light inhibiting amber bags that restrict exposure of the product to ultraviolet (UV) light.

Helapet Amber Bags are fully certified to Pharmacopoeia standards restricting the amount of light passing through the bag, not exceeding 10% threshold at any wavelength in the range of 290nm to 450nm.

Offering valuable UV protection for any healthcare application, Helapet Amber Bags are ideal for containing IV bags, pre-filled syringes and vials. Open-ended bag formats can be securely heat sealed for added security and resealable formats allow for multi-opening.

Helapet UV Light Inhibiting Amber Bags are supplied with supporting certification and lot traceability.



PB0192



Open Ended



Resealable



Sealable

Features and Benefits

- ✓ Strong 250 gauge / 62 micron polyethylene material
- ✓ Conforms to Pharmacopoeia standards*
- ✓ Multiple size options
- ✓ Open-ended, sealable and resealable
- ✓ Full lot traceability
- ✓ Non sterile

*All Helapet amber bags conform to US Pharmacopoeia performance testing standards for light transmission (661 & 671), not exceeding 10% threshold at any wavelength in the UV range of 290nm to 450nm.

Enquire today for **TRIAL** samples



Ordering Information

Open ended

PB0190

W203 x L355mm
non sterile • 500 pcs/pck

PB0191

W127 x L203mm
non sterile • 500 pcs/pck

PB0192

W76 x L228mm
non sterile • 1000 pcs/pck

PB0194

W152 x L355mm
non sterile • 1000 pcs/pck

PB0195

W279 x L457mm
non sterile • 500 pcs/pck

Re-sealable

PB0208

W152 x L203mm
non sterile • 1000 pcs/pck

PB0214

W203 x L355mm
non sterile • 1000 pcs/pck

Sealable

PB2012

W318 x L510mm + 36mm flap
non sterile • 500 pcs/pck

Technical Data

| Open ended | PB0190 | PB0191 | PB0192 | PB0194 | PB0195 |
|-----------------|---|-------------|--|-------------|---|
| Dimensions (mm) | W203 x L355 | W127 x L203 | W76 x L228 | W152 x L355 | W279 x L457 |
| Material | Virgin Low Density Polyethylene (LDPE) | | | | |
| Thickness | 250 gauge / 62 micron | | | | |
| Packaging | 100 x bags per resealable polybag; 5 x polybags per pack | | 100 x bags per resealable polybag; 10 x polybags per pack | | 100 x bags per resealable polybag; 5 x polybags per pack |
| Sterility | Non sterile | | | | |

Resealable

| | PB0208 | PB0214 |
|-----------------|--|-------------|
| Dimensions (mm) | W152 x L203 | W203 x L355 |
| Material | Virgin Low Density Polyethylene (LDPE) | |
| Thickness | 250 gauge / 62 micron | |
| Packaging | 100 x bags per resealable polybag; 10 x polybags per pack | |
| Sterility | Non sterile | |

Sealable

| | PB2012 |
|-----------------|---|
| Dimensions (mm) | W318 x L510mm + 36mm flap |
| Material | Virgin Low Density Polyethylene (LDPE) |
| Thickness | 200 gauge / 50 micron |
| Packaging | 100 x bags per resealable polybag; 5 x polybags per pack |
| Sterility | Non sterile |

Helapet Limited

Lister House | Blackburn Road | Houghton Regis | Bedfordshire | LU5 5BQ

Web: helapet.co.uk | Phone: +44 (0)1582 501980 | Email: sales@helapet.co.uk

Helapet Limited is a Clinimed Group Company.

