

UV Light Inhibiting Amber Bags



Certain medications, antibiotics and compounded chemotherapy products need to be stored and transported in light inhibiting amber bags that restrict exposure of the product to ultraviolet (UV) light.

Helapet Amber Bags are fully certified to Pharmacopoeia standards restricting the amount of light passing through the bag, not exceeding 10% threshold at any wavelength in the range of 290nm to 450nm.

Offering valuable UV protection for any healthcare application, Helapet Amber Bags are ideal for containing IV bags, pre-filled syringes and vials. Open-ended bag formats can be securely heat sealed for added security and resealable formats allow for multi-opening.

Helapet UV Light Inhibiting Amber Bags are supplied with supporting certification and lot traceability.



PB0192



Open Ended



Resealable



Sealable

Features and Benefits

- ✓ Strong 250 gauge / 62 micron polyethylene material
- ✓ Conforms to Pharmacopoeia standards*
- ✓ Multiple size options
- ✓ Open-ended, sealable and resealable
- ✓ Full lot traceability
- ✓ Non sterile

*All Helapet amber bags conform to US Pharmacopoeia performance testing standards for light transmission (661 & 671), not exceeding 10% threshold at any wavelength in the UV range of 290nm to 450nm.

Enquire today for
FREE
samples



Ordering Information

Open ended

PB0190

W203 x L355mm
non sterile • 500 pcs/pck

PB0191

W127 x L203mm
non sterile • 500 pcs/pck

PB0192

W76 x L228mm
non sterile • 1000 pcs/pck

PB0194

W152 x L355mm
non sterile • 1000 pcs/pck

PB0195

W279 x L457mm
non sterile • 500 pcs/pck

PB0196

W76 x L279mm
non sterile • 1000 pcs/pck

Re-sealable

PB0208

W152 x L203mm
non sterile • 1000 pcs/pck

PB0214

W203 x L355mm
non sterile • 1000 pcs/pck

Sealable

PB2012

W318 x L510mm + 36mm flap
non sterile • 500 pcs/pck

Technical Data

Open ended

	PB0190	PB0191	PB0192	PB0194	PB0195	PB0196
Dimensions (mm)	W203 x L355	W127 x L203	W76 x L228	W152 x L355	W279 x L457	W76 x L279
Material	Virgin Low Density Polyethylene (LDPE)					
Thickness	250 gauge / 62 micron					
Packaging	100 x bags per resealable polybag; 5 x polybags per pack		100 x bags per resealable polybag; 10 x polybags per pack		100 x bags per resealable polybag; 5 x polybags per pack	100 x bags per resealable polybag; 10 x polybags per pack
Sterility	Non sterile					

Resealable

	PB0208	PB0214
Dimensions (mm)	W152 x L203	W203 x L355
Material	Virgin Low Density Polyethylene (LDPE)	
Thickness	250 gauge / 62 micron	
Packaging	100 x bags per resealable polybag; 10 x polybags per pack	
Sterility	Non sterile	

Sealable

	PB2012
Dimensions (mm)	W318 x L510mm + 36mm flap
Material	Virgin Low Density Polyethylene (LDPE)
Thickness	200 gauge / 50 micron
Packaging	100 x bags per resealable polybag; 5 x polybags per pack
Sterility	Non sterile

