

# Vented Vial Adapter



Helapet Vented Vial Adapters are a safe and easy-to-use solution for reconstitution of cytotoxics, antibiotics and other hazardous chemicals.

A single-use needle-free device, Helapet Vented Vial Adapters can be securely fastened to any standard 13mm or 20mm diameter vial. The novel skirt design prevents the adapter from being detached once applied, providing irremovable and stable vial access. An integrated grip grants superior handling and reduced risk of contact with the luer lock connection.

Each device incorporates a bilateral 0.2µm hydrophobic filter allowing the user to precisely aspirate from the vial without the risk of aerosol contamination, or the need for manual pressure equalisation.

Helapet Vented Vial Adapters are presented individual blister packed, CE marked with full lot traceability



IV0013

IV0020L

## Features and Benefits

- ✓ Bilateral 0.2 µm vent filter
- ✓ Plastic spike for needle-free reconstitution
- ✓ Easy snap-on attachment
- ✓ Integrated grip reducing luer contact
- ✓ Compatible with cytotoxic substances
- ✓ Sterile
- ✓ CE marked

## Ordering Information

### IV0013

Vented Vial Adapter  
13mm • sterile • 100 pcs/pck

### IV0020L

Vented Vial Adapter  
20mm • sterile • 100 pcs/pck

Enquire today for  
**FREE**  
samples



## Technical Data

	IV0013	IV0020L
Housing material	MABS (Methyl Methacrylate Acrylonitrile Butadiene Styrene); LDPE (Low Density Polyethylene)	
Vent Filter material	Acrylic Copolymer / Nylon	
Adapter diameter	13mm	20mm
Air filter diameter	5mm	15mm
Attachment force	<27.7lbs	
Vent filter	0.2µm hydrophobic	
Flow rate	>40ml/min	>200ml/min
Priming volume	<0.08ml	
Luer access	Accepts ISO 594 Compliant Male Luer	
Packaging	Individually blister packed; 100 x blisters per pack	
Sterilisation	Gamma sterilised	
Compliance	CE marked; 510(k) cleared by United States Food and Drug Administration	
Shelf Life	3 years from date of manufacture	

Helapet Limited

Lister House | Blackburn Road | Houghton Regis | Bedfordshire | LU5 5BQ

**Web:** helapet.co.uk | **Phone:** +44 (0)1582 501980 | **Email:** sales@helapet.co.uk

Helapet Limited is a CliniMed Group Company.

