



Viscot Sterile Utility Marker Pen

The Viscot Sterile Utility Marker Pen has a durable fine tip, ideally suited for writing on medication labels within sterile working environments. Waterproof permanent ink supports writing on labels, plastic (i.e syringes) and glass without smudging or the risk of removal. Each pen includes 9 blank medication labels and a handy ruler, all supplied in sterile blister packaging.

Enquire today for **FREE** samples

Features and Benefits

- ✓ Permanent, non-smearing, non-toxic ink
- ✓ Ideal for labelling bags, vials and syringes
- ✓ Fine tip for precise writing
- ✓ Ink Colour: Black
- ✓ Includes 9 blank labels and ruler
- ✓ Sterile blister packed for cleanroom use



PT0141

Ordering Information

PT0141
sterile • 50 pcs/pck



Technical Data

Tip size	Fine
Ink Colour	Black
Properties	Permanent, non-smearing, non-toxic ink; latex free
Label size	L50mm x W13mm x 9 labels
Packaging	Individually blister packed; 50 x blisters per pack
Compliance	CE marked; Individually lot numbered
Sterilisation	Gamma sterilised
Shelf life	3 years from date of manufacture

Helapet Limited | Lister House | Blackburn Road | Houghton Regis | Bedfordshire | LU5 5BQ

Web: helapet.co.uk | **Phone:** +44 (0)1582 501980 | **Email:** sales@helapet.co.uk

Helapet Limited is a CliniMed Group Company.