



Scan Me for Web Page



PharmaPack[®]

sterile multi-pack syringes

Use PharmaPack[®] to safely and conveniently transfer several sterile syringes in a single bag, reducing time, energy, packaging waste and risk of contamination. Helapet offer a comprehensive range of quality prepared CE marked luer-lock and luer-slip syringes, available in several sizes and packaging configurations to suit specific production requirements.

Designed for use in aseptic manufacturing units, PharmaPack[®] supports the latest disinfection and aseptic transfer protocols with multi-bag packaging using low-lint materials. All packs are terminally sterilised for use in high grade cleanroom environments and undergo regular sample testing for added customer assurance.

PharmaPack[®] sterile BD luer lock syringes in combination with selected Helapet syringe caps are microbiologically and structurally integrity tested in accordance with the NHS Method Document: Protocols for the Integrity Testing of Devices*.



Features and Benefits

- ✓ Metric syringe graduation
- ✓ Optimised pack quantities
- ✓ Double or triple bagged
- ✓ Low-lint packaging materials
- ✓ Capped option to reduce contamination risk
- ✓ Sterilised for use in Grade A environments

Technical Data

Syringe Sizes	1ml, 3ml, 5ml, 10ml, 20ml, 30ml, 50ml	
Packaging Material	ETO sterilisation	Flexopeel [®] header pouches
	Gamma sterilisation	Clear polyethylene bags
Packaging Specifications	Sealed double and triple bagged formats	
Syringes per pack	5s, 10s, 25s, 50s	
Syringe Connections	Luer Lock; Luer Slip	
Capping	Pre capped with high density polyethylene white vented cap or blue non-vented cap (optional)	
Processing	ISO Class 7 rated cleanroom	
Compliance	Microbiologically and structurally integrity tested in accordance with the NHS Method Document: Protocols for the Integrity Testing of Devices*; CE-Marked; Bag seal integrity tested	
Sterilisation	Luer Lock 20ml, 30ml, 50ml	Sterilised by Ethylene Oxide (ETO)
	Luer Lock 1ml, 3ml, 5ml, 10ml + all triple bagged formats	Gamma irradiated to no less than 25 kGy
	Luer Slip 1ml, 2ml	Sterilised by Ethylene Oxide (ETO)
Shelf Life	3 years from date of manufacture	

Ordering Information

LUER LOCK / UNCAPPED

Double bagged packaging

XX3101BD

PharmaPack® BD syringe (LL)
1ml • Sterile • 50 x 10 pcs



XX3103BD

PharmaPack® BD syringe (LL)
3ml • Sterile • 50 x 10 pcs



XX3203BD

PharmaPack® BD syringe (LL)
3ml • Sterile • 20 x 25 pcs




XX3503BD

PharmaPack® BD syringe (LL)
3ml • Sterile • 10 x 50 pcs




XX3105BD

PharmaPack® BD syringe (LL)
5ml • Sterile • 40 x 10 pcs



XX3205BD

PharmaPack® BD syringe (LL)
5ml • Sterile • 16 x 25 pcs



XX3110BD

PharmaPack® BD syringe (LL)
10ml • Sterile • 35 x 10 pcs



XX3120BD

PharmaPack® BD syringe (LL)
20ml • Sterile • 10 x 10 pcs



XX3220BD

PharmaPack® BD syringe (LL)
20ml • Sterile • 4 x 25 pcs



XX3130BD

PharmaPack® BD syringe (LL)
30ml • Sterile • 10 x 10 pcs




XX3230BD

PharmaPack® BD syringe (LL)
30ml • Sterile • 4 x 25 pcs



XX3150BD

PharmaPack® BD syringe (LL)
50ml • Sterile • 5 x 10 pcs



XX3250BD

PharmaPack® BD syringe (LL)
50ml • Sterile • 2 x 25 pcs



XX3350BD

PharmaPack® BD syringe (LL)
50ml • Sterile • 5 x 25 pcs



Triple bagged packaging

XXTB30501BD

PharmaPack® BD syringe (LL)
1ml • Sterile • 50 x 5 pcs



XXTB30505BD

PharmaPack® BD syringe (LL)
5ml • Sterile • 50 x 5 pcs



XXTB30510BD

PharmaPack® BD syringe (LL)
10ml • Sterile • 40 x 5 pcs



XXTB3130BD

PharmaPack® BD syringe (LL)
30ml • Sterile • 8 x 10 pcs



XXTB3150BD

PharmaPack® BD syringe (LL)
50ml • Sterile • 4 x 10 pcs



Further double and triple bagged options available on request

LUER LOCK / CAPPED

Double bagged packaging


XX3101BDC

PharmaPack® BD syringe (LL)
1ml • Sterile • 50 x 10 pcs




XX3103BDC

PharmaPack® BD syringe (LL)
3ml • Sterile • 50 x 10 pcs




XX3203BDC

PharmaPack® BD syringe (LL)
3ml • Sterile • 20 x 25 pcs



XX3105BDC

PharmaPack® BD syringe (LL)
5ml • Sterile • 40 x 10 pcs



XX3205BDC

PharmaPack® BD syringe (LL)
5ml • Sterile • 16 x 25 pcs




XX3110BDC

PharmaPack® BD syringe (LL)
10ml • Sterile • 35 x 10 pcs



XX3120BDC

PharmaPack® BD syringe (LL)
20ml • Sterile • 10 x 10 pcs



XX3220BDC

PharmaPack® BD syringe (LL)
20ml • Sterile • 4 x 25 pcs



Triple bagged packaging

XXTB30505BDC

PharmaPack® BD syringe (LL)
5ml • Sterile • 50 x 5 pcs




XXTB30510BDC

PharmaPack® BD syringe (LL)
10ml • Sterile • 40 x 5 pcs



XXTB3120BDC

PharmaPack® BD syringe (LL)
20ml • Sterile • 8 x 10 pcs




LUER SLIP / UNCAPPED

Double bagged packaging

XX3101BBLS

PharmaPack® BB syringe (LS)
1ml • Sterile • 30 x 10 pcs



XX3101LS

PharmaPack® BD syringe (LS)
1ml • Sterile • 20 x 25 pcs



XX3052BD

PharmaPack® BD syringe (LS)
2ml • Sterile • 75 x 5 pcs



Helapet Syringe Caps

IV2102

Combi-Cap Luer Lock Syringe Cap
Sterile • 1,000 pcs/pck



IV5750R

Tamper Evident IV Syringe Cap
Sterile • 10 x 10 pcs/pck



PharmaPack® Integrated Blue Non-Vented Syringe Cap
(Non-stock item)



Syringe/endcap integrity tested in accordance with NHS Protocols for the Integrity Testing of Devices*

*NHS Method: Protocols for the Integrity Testing of Devices; 2nd Edition, revised April 2013

Helapet Limited

Lister House | Blackburn Road | Houghton Regis | Bedfordshire | LU5 5BQ

Web: helapet.co.uk | Phone: +44 (0)1582 501980 | Email: sales@helapet.co.uk

Helapet Limited is a Clinimed Group Company.

INVESTORS IN PEOPLE®
We invest in people Gold

