



Scan Me for Web Page



Sterile Single Use Goggles

Constructed from lightweight ultra-soft PVC with a toughened polycarbonate lens, BioClean™ Clearview Sterile Single Use Goggles are PPE Category 2 certified and designed to provide effective personal protection for users operating in controlled working environments. Non-linting and entirely latex-free, each pair offers assured protection against flying fragments, large chips, particles, chemical splashes and impacts. BioClean™ Clearview Sterile Single Use Goggles are optically correct with anti-fog and anti-scratch coating for clear vision. An indirect ventilation system protects against the risk of hairs and particles entering the cleanroom environment, as well as restricting potentially harmful material from reaching the interior of the goggles. Available in a universal size with adjustable head band, that can be comfortably worn over prescription eyewear.



Features and Benefits

- ✓ Lightweight
- ✓ Indirect ventilation system
- ✓ Anti-fog and scratch-resistant lens
- ✓ Non-linting and latex-free
- ✓ Individually packaged
- ✓ Sterilised by Ethylene Oxide (ETO)



CS0015

Ordering Information

CS0015

BioClean™ Sterile Single Use Goggles
sterile • 10 pcs/pck

Enquire
today for
FREE
samples



Technical Data

Material	Frame	Lightweight ultra-soft PVC
	Lens	Toughened polycarbonate
	Head Band	Latex free silicone
Material Weight	74g	
Lens Thickness	2mm	
Colour	Transparent frame and strap	
Properties	Non-linting, latex free construction	
Packaging	Individually packaged in sealed Tyvek™ / PE inner bag; 10 x goggles per pack	
Compatibility	ISO Class 4 and EU GMP Grade A	
Compliance	CE marked. Category 2 Personal Protective Equipment (PPE) Regulation (EU)2016/425. Complies with requirements of EN 166:2001 standard for personal eye protection	
Sterilisation	Ethylene Oxide (ETO)	
Shelf Life	3 years from date of manufacture Store in dry, cool place (<40°C) away from direct sunlight and fluorescent light	

BioClean™ is a registered trademark owned by Ansell Limited or one of its affiliates.

Helapet Limited

Lister House | Blackburn Road | Houghton Regis | Bedfordshire | LU5 5BQ

Web: helapet.co.uk | **Phone:** +44 (0)1582 501980 | **Email:** sales@helapet.co.uk

Helapet Limited is a CliniMed Group Company.

