



# berner

## ChemoSorb Pad

A highly absorbent pad which safely contains spilt liquids including cytotoxic spills. Liquids are converted into a gel which allows for safe and easy disposal.

### Features and Benefits

- ✓ Highly absorbent
- ✓ Non-drip binding of spilt liquids
- ✓ Suitable for cytotoxic drugs
- ✓ Low linting
- ✓ Latex free



### Ordering Information

**CS120016**  
 Berner ChemoSorb Pad  
 non sterile • 10 pcs/pck



### Technical Data

Material	Cellulose-fibre Two-way component polyolefin fibre
Colour	White
Dimensions (cm)	38cm x 54cm
Weight (kg)	0.105kg ± 7%
Packaging	10 x pads per pack
Absorbency (litres)	Up to 3 litres
Compatibility	Ideal for CMR (Carcinogenic, Mutagenic, toxic to Reproduction) drugs, for use in laboratories, cleanrooms and hospitals
Shelf life	3 years from the date of manufacture