

# Veterinary Adhesive Remover Spray



## Gentle, pain-free removal of wound dressings within seconds

Helapet Veterinary Adhesive Remover Spray is our ready-to-use solution for the quick, no-sting removal of wound dressings and medical devices using skin adhesives.

Repeated skin-stripping can progressively reduce the skin's integrity, especially sensitive areas that can result in unpleasant inflammation and irritation. This may lead animals to constantly bite and scratch resulting in more serious infections, impeding the healing process.

Helapet Veterinary Adhesive Remover spray grants comfortable, cold-free application to the skin reducing the risk of animal movement. Integrated "bag-on-valve" technology improves usability with quiet dispensing of the aerosol in any direction, preventing unnecessary stress to the animal during application.

Simple, pain-free removal is achieved by spraying around the edge of the adhesive and waiting just a few seconds before gently peeling back. The quick-drying formula leaves no residue that could compromise adhesion of replacement dressings, tape or medical devices, supporting quality care and reducing wound exposure.



## Features and Benefits

- ✓ Quick, no-sting removal of adhesive based dressings from skin, hair and fur
- ✓ Comfortable no-cold sensation when applied
- ✓ Reduced skin trauma and risk of infection resulting from animal bites and scratching
- ✓ Does not compromise adhesion once evaporated, allowing for precise repositioning or application of new dressings

## Our innovative "bag-on-valve" technology delivers

- ✓ Immediate use with no need to shake the can
- ✓ Quiet spray reducing risk of animal stress
- ✓ Consistent 360° spray distribution
- ✓ 100% dispensing of solution reducing wastage

## Ordering Information

**VAR1050**  
50ml • non sterile • 12 pcs/pck