



BD Plastipak™ 60ml Luer Lock Amber Syringe

Helapet offer the BD Plastipak™ 60ml Luer Lock Amber Syringe for the proven protection of photosensitive medications from production to administration.

Ideal for supporting compounding of photosensitive medications in aseptic and pharmacy manufacturing units, the BD Plastipak™ 60ml Luer Lock Amber Syringe has a distinctive coloured barrel to effectively block out up to 90% of light rays, keeping light sensitive drugs safe from breakdown or damage caused by UV, ambient and solar rays. Single 1ml graduated metric marking ensures precise dosages of medicines.

Presented individually sterile in blister packaging, offers the user safe and economical aseptic transfer into sterile working environments including cleanroom and isolators, reducing packaging waste and risk of cross contamination.



Features and Benefits

- ✓ Validated safe protection of photo sensitive drugs*
- ✓ 1ml metric graduation for precise dosage
- ✓ Secure fastening to luer lock devices and equipment
- ✓ Latex free construction reducing risk of allergic reaction
- ✓ Individually sterile blister packaged for economical use

*As per the transparency test (method 1) as described in the Japanese Pharmacopoeia XVI. According to such method the light transmissibility (%) is characterized under a UV light source emitting at 450nm

Ordering Information

BD3050AMB

BD Plastipak™ 60ml Luer Lock Amber Syringe
60ml • sterile • 4 x 60 pcs/pck

Enquire
today for
FREE
samples



Technical Data

Syringe Size	60ml												
Material	<table border="0"> <tr> <td>Barrels</td> <td>Polypropylene</td> </tr> <tr> <td>Plunger rods</td> <td>Polypropylene</td> </tr> <tr> <td>Barrels for BD catalogue number 309628</td> <td>Polycarbonate</td> </tr> <tr> <td>Stoppers</td> <td>Latex free elastomer</td> </tr> <tr> <td>Lubricant</td> <td>Medical grade silicone oil, <0.25mg/cm²</td> </tr> <tr> <td>Scale</td> <td>Ink and dissolvent</td> </tr> </table>	Barrels	Polypropylene	Plunger rods	Polypropylene	Barrels for BD catalogue number 309628	Polycarbonate	Stoppers	Latex free elastomer	Lubricant	Medical grade silicone oil, <0.25mg/cm ²	Scale	Ink and dissolvent
Barrels	Polypropylene												
Plunger rods	Polypropylene												
Barrels for BD catalogue number 309628	Polycarbonate												
Stoppers	Latex free elastomer												
Lubricant	Medical grade silicone oil, <0.25mg/cm ²												
Scale	Ink and dissolvent												
Colour	Opaque / Amber												
Packaging Material	Polyamide/polyethylene film; Medical grade paper												
Packaging Specifications	1 x individually sterile packaged syringe, 60 x syringes shrink wrapped per carton; 4 x cartons per case												
Syringe Connection	Luer Lock; Concentric												
Compliance	<p>CE certified Class IIa Medical Device under Rule 2, Annex IX of the Medical Devices Directive 93/42/EEC, AEMPS (0318) Certificate N°95 06 0005; ISO 13485 Certificate N° 2012 07 0013 EN</p> <p>Certified to reduce UV damage to light sensitive medications. The light transmission has been characterized as per the transparency test (method 1) as described in the Japanese Pharmacopeia XVI.</p> <p>According to such method the light transmissibility (%) is characterized under a UV light source emitting at 450nm. The light transmission is 5.6± 0.2% (mean ± standard deviation) according to the transparency test.</p>												
Sterilisation	Ethylene Oxide Sterilisation (ETO) following EN ISO 11135-1. EO residues are within applicable regulation												
Shelf Life	5 years. Store in a dry and warm place, not exposed to strong light												

Helapet Limited

Lister House | Blackburn Road | Houghton Regis | Bedfordshire | LU5 5BQ

Web: helapet.co.uk | **Phone:** +44 (0)1582 501980 | **Email:** sales@helapet.co.uk

Helapet Limited is a CliniMed Group Company.



CERTIFICATE NO. MD 78785