

Coloured Pharmacy Prep Trays

Sterile and Non Sterile

Helapet offer our range of single-use polypropylene prep trays available in four common colours, red, blue, yellow and purple for colour-coding applications*. Available sterile or non sterile in two convenient sizes, our trays are designed to support safe handling and transfer of consumables, instruments, medical devices and medications within cleanroom and isolator environments.

Scan Me for Web Page



Features and Benefits

- ✓ Lightweight, durable polypropylene material
- ✓ Four opaque colours for colour-coding applications and easy identification*
- ✓ Raised, smooth edges for safe containment and transportation
- ✓ Two convenient sizes for large and small consumable volumes
- ✓ Individually double or triple polybagged and terminally sterilised for assured material transfer



Technical Data

Tray Size	Large		Small	
Internal Base Dimensions (L x W x D)	295 x 245 x 55mm		180 x 115 x 35mm	
Material	Polypropylene			
Colours	Blue, Red, Yellow, Purple*			
Packaging	Sterile	Double Poly Bagged		
		Triple Poly Bagged		
Processing	Sterile trays are cleaned and packaged in an ISO Class 7 cleanroom			
Sterility	Sterile	Gamma irradiated to no less than 25 kGy		
	Non Sterile	Can be sterilised by gamma irradiation, ethylene oxide (ETO) and autoclave		
Shelf Life	Sterile	5 years from date of manufacture		

*other colours are available subject to minimum order quantity

For order information, please see overleaf.

Helapet Limited | Lister House | Blackburn Road | Houghton Regis | Bedfordshire | LU5 5BQ


Web: helapet.co.uk | **Phone:** +44 (0)1582 501980 | **Email:** sales@helapet.co.uk

Helapet Limited is a CliniMed Group Company.

Coloured Pharmacy Prep Trays

Sterile and Non Sterile

Ordering Information



Tray Size:
large

	Blue		Red		Yellow		Purple	
	Product Code	Unit Of Sales	Product Code	Unit Of Sales	Product Code	Unit Of Sales	Product Code	Unit Of Sales
Non Sterile	PP0312B	50 pcs/pck	PP0312R	50 pcs/pck	PP0312Y	50 pcs/pck	PP0312P	50 pcs/pck
Sterile Double-Bagged	XX0312B	25 pcs/pck	XX0312R	25 pcs/pck	XX0312Y	25 pcs/pck	XX0312P	25 pcs/pck
Sterile Triple-Bagged	XXTB0312B	20 pcs/pck	XXTB0312R	20 pcs/pck	XXTB0312Y	20 pcs/pck	XXTB0312P	20 pcs/pck



Tray Size:
small

	Blue		Red		Yellow		Purple	
	Product Code	Unit Of Sales	Product Code	Unit Of Sales	Product Code	Unit Of Sales	Product Code	Unit Of Sales
Non Sterile	PP0326B	200 pcs/pck	PP0326R	200 pcs/pck	PP0326Y	200 pcs/pck	PP0326P	200 pcs/pck
Sterile Double-Bagged	XX0326B	88 pcs/pck	XX0326R	88 pcs/pck	XX0326Y	88 pcs/pck	XX0326P	88 pcs/pck
Sterile Triple-Bagged	XXTB0326B	70 pcs/pck	XXTB0326R	70 pcs/pck	XXTB0326Y	70 pcs/pck	XXTB0326P	70 pcs/pck

*other colours are available subject to minimum order quantity

Helapet Limited

Lister House | Blackburn Road | Houghton Regis | Bedfordshire | LU5 5BQ

Web: helapet.co.uk | Phone: +44 (0)1582 501980 | Email: sales@helapet.co.uk

Helapet Limited is a CliniMed Group Company.



CERTIFICATE NO. MD 78785