



Tyvek® Overboots



Designed from low-linting, robust Tyvek® material each boot can be easily donned over standard footwear. Two integrated wrapping bands allow safe and assured fastening below the knee, maintaining freedom of movement.

Compatible up to ISO Class 4 & EU GMP Grade A operating conditions, Helapet Tyvek® Overboots are ideal for enhanced lower body protection against chemical splashes and particle contamination.

Features and Benefits

- ✓ Lace ties
- ✓ Slip retardant sole
- ✓ Low linting material
- ✓ Compatible up to ISO Class 4 environments

Ordering Information

AT0010

Tyvek® Overboots
std • non sterile • 100 pairs/pck



AT0010

Technical Data

Material	Tyvek® flashspun polyethylene
Colour	White
Thickness	40µm embossed
Material Weight	0.3kg
Packaging	100 x pairs polybagged
Compliance	Category 1 PPE
Sterilisation	Non sterile
Shelf Life	5 years from date of manufacture

