



Swabable Vial Adapter

Swabable Vial Access Devices grant safe, drip-free aspiration of solutions from vials, bottles and other containers. The non-vented Swabable Vial Adapter can access standard drug vials for needle-free drug preparation and administration. Using an integrated Arisure® Neutral Valve provides a microbiologically closed connection, preventing the risk of accidental leaks and spills. The novel shrouded spike design ensures easy attachment and assured vial security with 13mm and 20mm vials.

Features and Benefits

- ✓ Safe aspiration of drugs and other solutions, reducing the risk of contamination
- ✓ Safe, needlefree access to 13mm and 20mm standard vials, bottles and containers
- ✓ Microbiologically closed luer connection for assured product protection
- ✓ Low attachment force with novel skirt design facilitates simple and secure vial insertion
- ✓ Individually blister packed and sterile for use in aseptic environments



IV0033

IV0034

Ordering Information

IV0033

Swabable Vial Adapter
13mm • sterile • 100 pcs/pck

IV0034

Swabable Vial Adapter
20mm • sterile • 100 pcs/pck

Enquire
today for
TRIAL
samples



Technical Data

	IV0033	IV0034
Housing material	MABS (Methyl Methacrylate Acrylonitrile Butadiene Styrene); Polycarbonate, Silicone	
Adapter diameter	13mm	20mm
Properties	Disinfectable Needle-free Valve; Latex Free; DEHP Free	
Attachment force	<27.7lbf	
Flow rate	113.44ml/min	105.23ml/min
Residual volume	0.15ml	0.19ml
Luer access	Dry luer lock connection with Arisure® Neutral Valve	
Packaging	Individually blister packed; 100 x blisters per pack	
Sterilisation	Gamma sterilised	
Compliance	CE marked; 510(k) cleared by United States Food and Drug Administration	
Shelf Life	3 years from date of manufacture	

Helapet Limited

Lister House | Blackburn Road | Houghton Regis | Bedfordshire | LU5 5BQ

Web: helapet.co.uk | **Phone:** +44 (0)1582 501980 | **Email:** sales@helapet.co.uk

Helapet Limited is a CliniMed Group Company.

