



Nutricair™

Enteral feeding tubes

Adults



NCE506AU

Designed for the safe delivery of enteral nutrition, fluids and medication.

Nutricair™ Enteral Feeding Tubes support safer enteral feeding practices by helping reduce the risk of misconnections.

Benefits

- **Safer Connections**
ENFit® compliant to support safer enteral feeding and reduce the risk of misconnections
- **Clinical Flexibility**
MRI compatible
- **Position Verification**
Radiopaque design supports placement visualisation
- **Patient Comfort**
Atraumatic distal tip
- **Ease of Use**
Colour-coded sizing and depth markings



Scan for product information and sample enquiries.

Nutricair™ Enteral feeding tubes - Adults



NCE506AU

Features

- ENFit® male connector
- Transparent polyurethane tube
- Radiopaque strip containing 40% barium sulphate
- Closed distal atraumatic tip
- Four lateral eyes
- Numerical markings every 5 cm (20-100 cm)
- 110 cm tube length
- Colour-coded link
- Sterile, single-use device
- Individually blister packed

Regulatory & Compliance

- ENFit® connector system compliant with ISO 80369-3
- CE Marked Class IIa Medical Device (CE2797)
- Manufactured under an ISO 13485 certified quality management system
- Sterile single-use device
- MRI compatible
- Latex free
- Free from DEHP (DOP), BBP, DBP, DIBP and DHP
- REACH compliant

Technical Data

Tube Size Dimensions

CH/Fr	External Diameter mm	Internal Diameter mm
06	2.0	1.3
08	2.65	1.55
10	3.3	2.15
12	4.0	2.6
14	4.65	3.0
16	5.3	3.3
18	6.0	4.1
20	6.65	4.5

Ordering Details

Standard Feeding Tubes

Product Code	NHSSC NPC Code†	CH/Fr	Colour Code	Length cm	Pack Qty
NCE106AU	FAS85062	6	Light Green	110	20
NCE506AU	FAS85062	6	Light Green	130	20
NCE108AU	FAS85073	8	Blue	110	20
NCE508AU	FAS85073	8	Blue	130	20
NCE110AU	FAS85063	10	Dark Grey	110	20
NCE510AU	FAS85063	10	Dark Grey	130	20
NCE112AU	FAS85052	12	White	110	20
NCE512AU	FAS85052	12	White	130	20
NCE114AU	FAS85070	14	Dark Green	110	20
NCE514AU	FAS85070	14	Dark Green	130	20
NCE116AU	FAS85056	16	Orange	110	20
NCE516AU	FAS85056	16	Orange	130	20
NCE118AU	-	18	Red	110	20
NCE518AU	-	18	Red	130	20
NCE120AU	-	20	Yellow	110	20
NCE520AU	-	20	Yellow	130	20

† National Product Code

Enteral Feeding Tubes with Guide-Wire

Product Code	NHSSC NPC Code†	CH/Fr	Colour Code	Length cm	Pack Qty
NCE106G	FAS85054	6	Light Green	110	20
NCE506G	FAS85054	6	Light Green	130	20
NCE108G	FAS85060	8	Blue	110	20
NCE508G	FAS85060	8	Blue	130	20
NCE110G	FAS85066	10	Dark Grey	110	20
NCE510G	FAS85066	10	Dark Grey	130	20
NCE112G	FAS85061	12	White	110	20
NCE512G	FAS85061	12	White	130	20
NCE114G	FAS85074	14	Dark Green	110	20
NCE514G	FAS85074	14	Dark Green	130	20
NCE116G	FAS85072	16	Orange	110	20
NCE516G	FAS85072	16	Orange	130	20

† National Product Code

Helapet Limited.

Lister House, Blackburn Road, Houghton Regis, Bedfordshire LU5 5BQ.

Phone: +44 (0)1582 501980 Email: sales@helapet.co.uk

helapet.co.uk



CERTIFICATE NO. MD 78785

Helapet Limited is a CliniMed Group Company. SteriClean® is a registered trademark of CliniMed (Holdings) Limited.