



# Arisure® Closed Male Luer

The Arisure® Closed Male Luer (CML) is easy to use and completely needle-free. It forms a dry disconnect with any Arisure® needle-free valve. It remains closed when not accessed, preventing accidental spills or discharge of hazardous drugs.

The Arisure® Closed Male Luer forms part of the Arisure® Closed System Transfer System, designed for healthcare professionals to ensure the safe compounding, transport, and delivery of hazardous drugs. This keeps you focused on what matters — the patient.

## Features and Benefits

- ✓ Dry disconnect with Arisure® needle-free valve for reduced exposure to harmful drug residue
- ✓ Safe and simple transfer of compounds
- ✓ Mechanically closed when not accessed, preventing accidental leaks or spills
- ✓ Microbiologically closed, reducing the risk of contamination
- ✓ Maintains seal integrity, keeping all connections safe even when accessed multiple times
- ✓ No special adapters means less time searching for parts and more time caring for your patients



IV0058

## Ordering Information

### IV0058

Arisure® Closed Male Luer  
sterile • 100 pcs/pck

Enquire  
today for  
**TRIAL**  
samples



## Technical Data

	IV0058
Housing material	Polycarbonate, Silicone
Vent filter material	Acrylic Copolymer / Nylon
Length (mm)	28.4
Properties	Disinfectable Needle-free Valve; Latex Free; DEHP Free
Air filter diameter	5mm
Attachment force	<27.7lbs
Vent filter	0.2µm hydrophobic
Flow rate	≥ 66 mL/minn
Residual volume	~ 0.05 mL
Luer access	Dry luer lock connection with Arisure® Neutral Valve
Packaging	
Sterilisation	Gamma sterilised
Compliance	CE marked; 510(k) cleared by United States Food and Drug Administration
Shelf Life	3 years from date of manufacture

## The Complete Arisure Closed Transfer System



IV0051 / IV0053  
Arisure® Closed Vial Adapters



IV0058  
Arisure® Closed Male Luer



IV0050  
Arisure® Dry Spike

### Helapet Limited

Lister House | Blackburn Road | Houghton Regis | Bedfordshire | LU5 5BQ

**Web:** helapet.co.uk | **Phone:** +44 (0)1582 501980 | **Email:** sales@helapet.co.uk

Helapet Limited is a CliniMed Group Company.

