



Arisure[®] Closed Vial Adapters

The Arisure[®] CVA is made to protect you from hazardous drug exposure. It is completely closed, mechanically trapping anything from escaping the vial and preventing contaminants from entering the vial. It is a safe and simple solution for reconstitution and transfer of drugs.

Available in both 13mm and 20mm sizes, the Arisure[®] Closed Vial Adapter forms part of the Arisure[®] Closed System Transfer System, designed for healthcare professionals to ensure the safe compounding, transport, and delivery of hazardous drugs. This keeps you focused on what matters — the patient.

Features and Benefits

- ✓ Mechanically closed to prevent the escape of harmful vapor and aerosols
- ✓ Reduced exposure to harmful drug residue
- ✓ Safe and easy transfer of compounds
- ✓ Microbiologically closed maintaining drug vial sterility for up to 7 days*
- ✓ Low attachment force and non-coring spike help eliminate problems during vial insertion
- ✓ CE Marked and Individually blister packed

*data available on request



IV0051 - 13mm

IV0053 - 20mm

Ordering Information

IV0051

Arisure[®] Closed Vial Adapter
13mm • sterile • 50 pcs/pck

IV0053

Arisure[®] Closed Vial Adapter
20mm • sterile • 50 pcs/pck

Enquire
today for
TRIAL
samples



Technical Data

	IV0051	IV0053
Material	Methyl Methacrylate; Acrylonitrile Butadiene Styrene; Polycarbonate; Thermo[plastic; Silicone	
Vent filter material	Polytetrafluoroethylene	
Adapter diameter	13mm	20mm
Drug Compatibility	Performance evaluated with a list of common hazardous medications*	
Air filter diameter	5mm	15mm
Attachment force	<20lbs	<40lbs
Vent filter	0.2µm hydrophobic and oleophobic	
Flow rate	>105ml/min	>108ml/min
Residual volume	0.2ml	0.3ml
Luer access	Compatible with male luers; Arisure® Closed Male Luer required for dry disconnect	
Packaging	Individually blister packed; 50 x blisters per pack	
Sterilisation	Gamma sterilised	
Compliance	CE marked; 510(k) cleared by United States Food and Drug Administration	
Shelf Life	3 years from date of manufacture	

*Data available on request

The Complete Arisure® Closed Transfer System



IV0051 / IV0053
Arisure® Closed Vial Adapters



IV0058
Arisure® Closed Male Luer



IV0050
Arisure® Dry Spike

Helapet Limited

Lister House | Blackburn Road | Houghton Regis | Bedfordshire | LU5 5BQ

Web: helapet.co.uk | **Phone:** +44 (0)1582 501980 | **Email:** sales@helapet.co.uk

Helapet Limited is a CliniMed Group Company.

